Product Information: AgriStrip

General Information

AgriStrip - the rapid one-step assay for the detection of plant pathogens

Intended use of the AgriStrip

The **AgriStrip** is intended to confirm the presence of a specific pathogen on-site in suspicious plant samples within minutes.

Test principle

AgriStrip, developed and produced by BIOREBA, is based on lateral flow immunochromatography using specific antibodies. Once the strip is dipped with the «sample» side into the sample extract, the liquid migrates upwards and initiates the antigen-antibody reaction which results in visible lines. Both test and control lines become visible with positive extracts, whereas negative samples produce the upper control line only (Fig. 5, next page). Lines start developing after 1-2 minutes and reach maximum intensity after 10-15 minutes. Dried test strips can be kept as permanent records.

Assay procedure:

- 1) Take/cut a part of a symptomatic plant sample and place it into an extraction bag (Fig. 1). Add extraction buffer with a disposable pipette. Please refer to the specific Product information for detailed instructions on **sample size, buffer type and sample/buffer ratio**.
- 2) Homogenize the sample with a handheld homogenizer on a flat surface (Fig. 2)
- 3) Transfer a total of 4 drops (approx. 150 µl) of the extract into a disposable cuvette (Fig. 3). Please refer to the specific Product information for **recommended dilutions of the extract**.
- 4) Insert one strip with the end marked «sample» into the extract and observe formation of colored bands (Fig. 4).

Fig. 1

Fig. 2

Fig. 3


Fig. 4


Interpretation of results

- 1) Both a strong test line and a strong control line (as shown in Fig. 5) are obtained with samples containing a high pathogen concentration.
- 2) A test line which appears weaker than the control line is typical for samples with low antigen (pathogen) concentration.
- 3) If only the control line appears, no detectable antigen is present in the sample. A very faint test line occasionally observed after 10-15 minutes should be interpreted as negative reaction.
- 4) If neither test line nor control line become visible, the test is invalid and should be repeated with a fresh strip.
- 5) A green test line without red component, which might be occasionally observed, should be interpreted as negative result. Repeat the test with a more diluted sample.

Note:

- 1) Store the strips and the extraction buffer at 4°C. However, exposure to ambient temperature (10 - 30°C), such as during transport and use in the field, does not affect the quality. Keep the packaging (containing desiccant bags) always hermetically closed. Absorbed moisture by the strips can lead to poor results or even complete failure of the test.
- 2) Strips must be used before the expiration date indicated on the label of the packaging.
- 3) Not all AgriStrip are run with the same extraction buffer or the same dilution. Refer to the indications on the specific Product information and use the buffer which is part of the complete kit.

Fig. 5: Test strips

Ordering information

The AgriStrip is available in three different ordering options: as complete kit (Fig. 6) or in bulk with 25 or 100 strips, respectively. Please refer to the internet www.bioreba.com for an updated list of tests available in the AgriStrip format.

Fig. 6: AgriStrip Complete kit (25 assays)

Content of the Complete kit

- ☒ 25 strips packed in a tube with desiccant bags
- ☒ 100 ml AgriStrip Extraction buffer, ready-to-use
- ☒ 25 disposable pipettes
- ☒ 25 extraction bags «Universal»
- ☒ 25 disposable cuvettes
- ☒ 1 cuvette rack

Note: the hand-held homogenizer (Art. No. 40010) is optional and not included in the complete kit.

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