Product Information: AgriStrip-magnetic Grapevine leafroll-assoc. virus 3 (GLRaV-3)

GLRaV-3 AgriStrip-magnetic - a rapid assay for the detection of Grapevine leafroll-associated virus 3 (GLRaV-3)

Intended use

The rapid assay GLRaV-3 AgriStrip-magnetic from BIOREBA enables the identification of GLRaV-3 in symptomatic tissue of grapevine. Optimally, the presence of GLRaV-3 is assessed in mature grapevine leaves in autumn or in bark scrapings during dormancy. It is not recommended to use young tissue early in the growing season.

The symptoms of GLRaV-3 cannot be distinguished from those caused by other grapevine viruses. Red cultivars show a reddish leaf area, whereas major veins stay green. White cultivars typically show slight yellowish discoloration of leaves. Furthermore, the leaves roll characteristically.

Test principle

The innovative rapid assay GLRaV-3 AgriStrip-magnetic, developed and produced by BIOREBA, is a lateral-flow immunochromatographic test based on an antigen-antibody reaction combined with Magnetic Bead Technology, enabling the preceding enrichment of the antigen. Grapevine leaf extract or extract of bark scrapings is first incubated with magnetic beads coated with antibodies specific for GLRaV-3. These beads are then separated from the extract with a magnet and resuspended in a small volume of running buffer. Finally, the concentrated beads migrate upwards through the test strip and brown colored lines become visible.

Both test and control lines become visible with positive extracts (containing GLRaV-3), whereas negative samples produce only the upper control line (Fig. A). Intense coloration is reached within 10-20 min. and the result can be registered.


Specificity and sensitivity

The GLRaV-3 AgriStrip-magnetic is based on monoclonal antibodies and reacts specifically with GLRaV-3. No cross-reaction is observed with other grapevine viruses. Strongest coloration of test bands is obtained with a 1:30 (w/v) dilution of an extract of infected grapevine leaves or bark scrapings.

The sensitivity attained with the AgriStrip lateral flow is comparable to the sensitivity in the DAS-ELISA method. Evaluation studies of several laboratories involved in grapevine certification revealed an accordance of 100% between the two test formats with field samples.

Please read the notes on page 3 before starting the assay.

Fig. A
Test strip

Test procedure

Sampling and extraction


Label extraction bags, tubes and strips with a permanent marker to prevent confusion.

1) Place a piece of **0.1 g** of symptomatic leaf (that is approximately the size of a 1-EURO coin) or bark scrapings in an extraction bag and add 3 ml of **AgriStrip Extraction buffer C** with a disposable pipette*.

2) Homogenize the tissue with the homogenizer hand model on a flat surface with a few movements for 5-7 sec (crushing bark scrapings first with a hammer is recommended) to obtain a 1:30 (w/v) extract.

3) Remove the lid of a GLRaV-3 magnetic beads tube. Place the tube at a **front position** (not magnetic) of the rack. Insert the tube as deep as possible. With a disposable pipette*,

transfer 1 ml of the extract from the extraction bag into the tube (fill it until the liquid layer becomes visible in the rack). Close the tube.

4) Remove the tube from the rack and shake it (by turning upside down) until the brownish content is completely resuspended. Place the tube again at a **front position** (not magnetic) on the rack and wait for **7 min** ⌚

5) Remove the tube and shake it for 3 sec. Remove the lid. Insert the tube as deep as possible at the position «**magnetic**» of the rack.

Wait for **3 min** ⌚ to allow the beads to accumulate at the wall of the tube by the effect of the magnet. Do not remove or turn the tube during this period of time in order not to disturb this process.


6) Discard the liquid in a waste container.

Addition of buffer

7) Remove the tube from the position «magnetic». A brownish precipitate is visible at the wall of the tube. Add **4 drops of Running buffer R** and close the lid.

Resuspension

8) **Shake vigorously for 10 sec** to resuspend the accumulated beads to obtain a homogenous suspension.

Immuno-chromatography

9) Place the tube at a **front position** (not magnetic) on the rack and insert a test-strip with the end marked «sample». Observe the coloration of the lines that is completed after 10-20 min.

Analysis of results

- 1) Both, a clearly visible test and control line are obtained with extracts containing GLRaV-3 as shown in Fig. A. This means the test is **positive**.
- 2) The test is **negative**, if the control line but no test line appears after 10-20 min. This means that

no GLRaV-3 is in the extract or in a concentration below the detection limit.

- 3) If neither test nor control line become visible, the test is **invalid** and should be repeated with a fresh strip.
- 4) Rarely, a faint test line might become visible after 20-30 min. Repeat the test in this case with a new extract or send the sample to a laboratory for verification with another method.

Dried test strips can be kept as permanent records, even though the coloration of bands is stronger when the strips are still wet.

Notes

- 1) Be aware that the strong magnets of the magnetic rack may interfere with or even destroy magnetic sensitive devices such as credit cards or watches.
- 2) Store the strips and the buffers at 4°C. However, limited exposure to ambient temperature (10 - 30°C), such as during transport and use in the field, does not affect the quality. Keep the packaging of the strips (containing desiccant bags) and the tubes containing the magnetic beads always hermetically closed. Absorbed moisture can lead to poor results or even complete failure of the test.
- 3) Strips must be used before the expiration date indicated on the label of the packaging.
- 4) Use the GLRaV-3 AgriStrip-magnetic always with the Extraction buffer C and the Running buffer R supplied with the strips. Other buffers might fail.
- 5) It is recommended to run the test at ambient temperature (15–25°C). At temperatures below 15°C, the speed of color development of test and control lines slows down.
- 6) The extraction buffer and running buffer contain sodium azide as preservative. Keep out of reach of children and do not ingest.

* Use of the disposable pipettes

Do not accidentally contaminate the extraction buffer. Please make sure that the disposable pipettes are immediately disposed of, after having been in contact with samples. Use always a **fresh pipette** from the pouch for dosing extraction buffer.

GLRaV-3 AgriStrip-magnetic

Fig. GLRaV-3 AgriStrip-magnetic Complete kit (Art. No. 122381)

Content of the Complete kit

- 25 strips packed with desiccant bags
- 25 tubes containing magnetic beads coupled to GLRaV-3 specific antibody
- 100 ml AgriStrip Extraction buffer C
- 5 ml Running buffer R
- 25 disposable pipettes
- 25 extraction bags
- Note: the magnetic rack and the homogenizer hand model are optional and not included in the Complete kit.

Ordering Information

Product	Art. No.	Size	Assays
GLRaV-3 AgriStrip-magnetic Complete kit	122381	1 Kit	25
<i>Optional tools required for the test but not supplied in the Complete kit</i>			
Magnetic rack	2367	1 unit	
Homogenizer hand model	400010	1 unit	

Disclaimer

This BIOREBA product is guaranteed to meet the specifications mentioned in this product information. No further warranties are given. Should this product fail for reasons other than inappropriate handling or misuse, BIOREBA AG will replace the product free of charge or refund the purchase price.

BIOREBA AG shall not be liable for any direct or indirect, special or consequential damages of any kind resulting from the use of this product.

Produktinformation: AgriStrip-magnetic Grapevine leafroll-assoc. virus 3 (GLRaV-3)

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