## Technical Information

# «Caulimo» buffer preparation

**Extraction buffer «Caulimo»:** 250 ml or 500 ml concentrate (10x) packed with powder of urea and egg albumin as package provided with the complete kits. Available also separately (Art. No. 110124).

### Preparation:

Make up 100 ml of the 10x concentrate to 1000 ml with distilled water. Then add slowly 60 g of urea while stirring until completely dissolved. This buffer is stable for 3 months at 4°C.

Just before use, add 10 g per litre egg albumin. First dissolve 10 g of egg albumin in about 20 ml of buffer and then make up to a total of 1000 ml with the remaining buffer (or at equal ratio for other volumes).

### Composition of the «Caulimo» extraction buffer:

for 1000 ml, in dist. water (pH 7.4 adjusted with HCl)

TRIS	2.4 g
NaCl	8.0 g
KCl	0.2 g
Tween 20	0.5 g
Triton X-100	5.0 g
PVP	10.0 g
Urea	60.0 g
NaN <sub>3</sub>	0.2 g
BSA	10.0 g
*Egg albumin (ovalbumin)	10.0 g

\* to be added just before use

### Biohazard

To prevent microbial growth, the buffers contain 0.02 - 0.2 % sodium azide which is toxic when swallowed and irritant if it comes in contact with skin, eyes and mucosae; nevertheless, these products are not considered to present a hazard to health when used in accordance to the principles of Good Laboratory Practice and Good Standards of Occupational Hygiene. Used reagent can be disposed of in the waste water in accordance with local regulations.

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