Product Information

Positive and negative controls

General description

Positive and negative controls are intended to verify the ELISA performance. The lyophilized positive controls are made from extracts of infected plants (virus reagents), or from extracts of heat-treated bacterial or mycelial suspensions (bacteria and fungi reagents). The negative controls represent extracts of healthy plants. The dry weight of the packed control is approx. 50-100 mg.

The positive controls are optimized for use in ELISA and are not suitable for other applications, e.g. inoculation of plants or our lateral flow test AgriStrip.

Note: Positive controls may contain more than one pathogen, hence a control may also give a positive reaction with more than one reagent.

Storage and use

Store the lyophilized controls at +4°C. The antigenic activity may gradually decrease somewhat over time, but the control reagent is usable until the expiration date marked on the datasheet.

Before use, add 2.5 ml distilled water (if not stated otherwise*) to the content and mix well; use this reconstituted extract immediately. The unused portion can be divided into aliquots and frozen at -20°C or colder for later use. However, this may result in a weaker response. Do not thaw until just before use and avoid refreezing.

* For testing potygroup (PTA-ELISA, Art. No. 163072/75), use the «PTA Extraction Buffer» instead of distilled water for reconstitution of controls.

Biosafety

Potential remaining infectivity is usually lost quickly after reconstitution of the lyophilized control. Even if the antigenic structure is still preserved to serve as controls in ELISA, there is realistically no danger for any transmission to crops if the controls are properly handled in laboratory tests as they are intended to.

Exp. date

The quality of BIOREBA positive and negative controls is always tested by ELISA. The expected values for the controls are indicated in the corresponding ELISA datasheet.

Positive and negative controls are stable for at least 12 months and up to 24 months from the delivery date. During this period, we guarantee the functionality of the controls in ELISA tests (carried out according to the BIOREBA procedure).

The validity of the positive and negative controls indicated on the datasheet is not specific to the lot number, unlike conjugates and IgGs. Therefore it is possible that positive or negative controls have the same lot number with a different expiration date.

Our negative control is used universally for all ELISA tests. If the negative control is ordered separately, no datasheet will be provided. In this case, BIOREBA guarantees a validity of 12 months from the date of delivery.

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