## AgriStrip<sup>cassette</sup> – the rapid one-step assay for the detection of plant pathogens

### Intended use


The AgriStrip<sup>cassette</sup> is intended to confirm the presence of a specific pathogen on-site in suspicious plant samples within minutes.

### Test principle

AgriStrip<sup>cassette</sup>, developed and produced by BIOREBA, is based on lateral flow immunochromatography using specific antibodies. Once a few drops of the sample extract are squeezed onto the sample spot of the test cassette, the liquid containing the antigens migrates initiating the antigen-antibody reaction which results in visible lines. Both test and control lines become visible with positive extracts, whereas negative samples produce the upper control line only (Fig. 4). Lines start developing after 1-2 minutes and reach maximum intensity after 15-20 minutes.

### Assay procedure

1. Remove the dropper tip from the provided homogenization tube. Take/cut a part of a symptomatic plant sample and place it into the tube containing the specific amount of extraction buffer and glass beads for homogenization (Fig. 1). Please refer to the specific Product Information for detailed instructions on sample size, buffer type, and sample/buffer ratio.
2. Homogenize the sample by shaking the tube thoroughly by hand and re-insert the dropper tip in the tube (Fig. 2).
3. Transfer 4 drops of the extract onto the sample spot of the test cassette (Fig. 3).
4. Observe the formation of test/control lines (Fig. 4).

## Interpretation of results

1. Both a strong test line and a strong control line (as shown in Fig. 4) are obtained with samples containing a high pathogen concentration.
2. A test line that appears weaker than the control line is typical for samples with low antigen (pathogen) concentration.
3. If only the control line appears (as shown in Fig. 4), no detectable antigen is present in the sample. A very faint test line occasionally observed after 10-15 minutes should be interpreted as a negative reaction.
4. If neither test line nor control line becomes visible, the test is invalid and should be repeated with a fresh strip.
5. A green test line without a red component, which might be occasionally observed, should be interpreted as a negative result. Repeat the test with a more diluted sample.

Fig. 5

### Note:

1. Store the unused cassettes and tubes containing extraction buffer at 4°C. Exposure to ambient temperature (10 - 30°C), such as during transport and use in the field, does not affect the quality of the product. Keep the packaging (containing desiccant bags) always hermetically closed. Absorbed moisture by the strips can lead to poor results or even complete failure of the test.
2. Cassettes and tubes must be used before the expiration date indicated on the label of the packaging.
3. Not all AgriStrip<sup>cassette</sup> are run with the same extraction buffer. Refer to the indications on the specific Product Information.

## Ordering Information

The AgriStrip<sup>cassette</sup> is available as a kit only, consisting of 5 assays (Fig. 5).

### Content of AgriStrip<sup>cassette</sup> kit

- ✓ 5 cassettes packed with desiccant bags
- ✓ 5 tubes containing glass beads and extraction buffer
- ✓ 5 dropper tips
- ✓ 1 technical information
- ✓ 1 product information

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