# Product Information

## Complete kits

### Kit content

| Component                   | 96 assays        | 480 assays        | 960 assays         |
|-----------------------------|------------------|-------------------|--------------------|
| IgG                         | 0.025 ml         | 0.1 ml            | 0.2 ml             |
| Conjugate                   | 0.025 ml         | 0.1 ml            | 0.2 ml             |
| Positive control            | 2.5 ml           | 2.5 ml            | 2.5 ml             |
| Negative control            | 2.5 ml           | 2.5 ml            | 2.5 ml             |
| Extraction buffer (10x/5x)* | 50 ml/2x50 ml    | 250 ml/500 ml     | 500 ml/2x500 ml    |
| Coating buffer              | 1 tablet         | 1 tablet          | 2 tablets          |
| Conjugate buffer (10x)      | 2 ml             | 10 ml             | 20 ml              |
| Substrate buffer (5x)       | 4 ml             | 20 ml             | 2 x 20 ml          |
| Washing buffer              | 1 tablet         | 1 pouch           | 2 pouches          |
| Substrate (pNPP)            | 4 tablets (5 mg) | 5 tablets (20 mg) | 10 tablets (20 mg) |
| Egg albumin**               | 5 g              | 25 g              | 50 g               |
| Urea***                     | 30 g             | 150 g             | 300 g              |
| Microtiter plates           | 1 x F96 plates   | 5 x F96 plates    | 10 x F96 plates    |

\* For grapevine testing, twice as much buffer is included since the concentration is 5x.

\*\* Is included only with the extraction buffer «Bulbs & Tubers» and «Caulimo».

\*\*\* Is included only with the extraction buffer «Caulimo».

### Intended use

All components of this kit are intended for use in combination with ELISA reagents for the detection of plant pathogens and not for use in human or animal application. For laboratory use only.

### Storage and stability

All components of the kit are stable until the indicated expiration date when stored at 4 °C. Once diluted, the buffers are stable for up to one month at 4 °C. Exceptions are the washing buffer which should be used within 2 days, and the extraction buffer containing egg albumin which has to be prepared freshly before use. Use all buffers at room temperature, i.e. it is advisable to bring the aliquot needed for the actual test to room temperature before use. **Note:** For long-term storage (>3 months) keep the substrate tablets at -20 °C.

1/1

Version: 6 - 17.03.2025

Adaption from last version: Update of BIOREBA and IQ Net logo.

[www.bioreba.com](http://www.bioreba.com)

**BIOREBA AG**

Christoph Merian-Ring 7  
CH-4153 Reinach  
Switzerland

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