Technical Information

Double antibody sandwich enzyme-linked immunosorbent assay (DAS-ELISA): test specifications

BIOREBA ELISA reagents were developed and optimized for application in the DAS-ELISA format (double antibody sandwich enzyme-linked immunosorbent assay).

General information

- All available ELISA products are optimized and validated by using certified Nunc MaxiSorp F96 microtiter plates.
- Optimal results are obtained by using certified Nunc MaxiSorp microtiter plates with a working volume of 200 µl per well.
- Buffers should be brought to room temperature before use.
- Laboratory vessels must be exclusively polyethylene (PE) or glass.
- After pipetting, the IgG and conjugate tubes must be stored immediately at 4 °C.
- Antibodies undergo denaturation when the wells are buffer-free between steps. After each washing the wells must be refilled within 10 minutes.
- BIOREBA recommends including a positive and negative control with each test.
- If any of the wells are not used in the test, please ensure that the empty wells do not come into contact with washing buffer. Wells that have been filled with washing buffer cannot be used for further tests.
- Avoid contact with skin when working with the ELISA kit components.
- Please note: The ELISA test for the «Poty group» is in PTA-format (Plate-trapped Antigen).

Recommendations for samples preparation

- For information on how to prepare the samples (sample dilution / appropriate extraction buffer, etc.) refer to the product information (PI) that is included with every delivery. The PI can also be downloaded as a pdf document at www.bioreba.com.
- The HOMEX homogenizer in combination with BIOREBA extraction bags is recommended for sample preparation.

Recommendations for the washing steps

- Initiate the washing only when the subsequent step in the procedure has been prepared.
- If it is not possible to refill the wells right after washing, the plate can be inverted immediately onto a wet, clean paper towel. After the «coating step» and after washing and taping the plate dry, it can be stored in a closed plastic bag for several weeks at -20 °C (does not apply to other steps in the process). Avoid repeated freeze thawing.

Alternative incubation conditions

Our optimized DAS-ELISA protocol is described on the next page. Nevertheless, there are alternative incubation conditions that will give similar results.

Coating step: A shorter incubation time such as 2.5 h at 37 °C gives similar results. Alternatively, an overnight incubation in the refrigerator at 4 °C is possible. Coated plates can be kept in the refrigerator at 4 °C for a maximum of one week (plates covered, preservative such as NaN₃ in the buffer).

Conjugation step: A shorter incubation time at 37 °C for 2 h gives similar results if the substrate incubation is prolonged to 1.5 h as compared to 1 h as standard (see next page).

Procedure


1. Coating: binding specific antibodies

- 1.1. Dilute the IgG 1:1000 in coating buffer (e.g. 20 µl IgG in 20 ml coating buffer).
- 1.2. Pipette 200 µl diluted IgG to each well.
- 1.3. Cover the plate (e.g. Parafilm or plate cover).
- 1.4. Incubate the plate at 30 °C for 4 h or at 4 °C overnight.
- 1.5. Wash the plate 3 x with 200 µl washing buffer per well (see point 5. and recommendations on the front page).


2. Antigen: incubation with plant extract

- 2.1. Homogenize the samples in extraction buffer (see recommendations on the front page).
- 2.2. Add 200 µl of the plant extract to each well.
- 2.3. Cover the plates with Parafilm or a plate cover and incubate at 4 °C overnight.
- 2.4. Wash the plate 3 x with 200 µl washing buffer per well (see point 5. and recommendations on the front page).

3. Conjugate: incubate with an enzyme-labelled antibody

- 3.1. Dilute the conjugate 1:1000 in conjugate buffer (e.g., 20 µl conjugate in 20 ml conjugate buffer).
- 3.2. Add 200 µl to each well.
- 3.3. Cover the plate with Parafilm or a plate cover and incubate at 30 °C for 5 h.
- 3.4. Wash the plate 3 x with 200 µl washing buffer per well (see point 5. and recommendations on the front page).

4. Substrate: the colour reaction indicates positive samples

- 4.1. Dissolve the substrate pNPP (para-nitrophenylphosphate) in substrate buffer to give a final concentration of 1 mg/ml.
- 4.2. Add 200 µl of the substrate solution to each well.
- 4.3. Incubate the plate at room temperature (20–25 °C in the dark).
- 4.4. Monitor the colour development visually and / or photometrically at 405 nm (for ELISA Readers with individual filters) or at 405/492 nm (for ELISA Readers with dual filters). Measure after 30–120 min.

5. Washing steps (washing the plate)

- 5.1. For the washing procedure empty the content of the plate by taping it dry. Ensure that there is nothing left in the well to avoid contamination from one well to the other.
- 5.2. Wash each well 3 times with 200 µl washing buffer. After each washing empty the contents of the plate by taping it dry.
- 5.3. After the last washing remove the remaining liquid by taping it onto a clean paper towel.

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