Product Information: AgriStrip

Spongospora subterranea (Sss)

Sss AgriStrip - a rapid assay for the detection of *Spongospora subterranea* f. sp. *subterranea* (Sss), and diagnosis of powdery scab on potato tuber

Intended use

The rapid assay Sss AgriStrip for *Spongospora subterranea* f. sp. *subterranea* detection is produced by BIOREBA and has been validated in cooperation with the Swiss Federal Institute of Technology (SFIT/ETH), and the Swiss College of Agriculture (SHL). This project was supported by the Innovation Promotion Agency (CTI) of the Swiss government.

The rapid assay enables confirmation of the presence of the powdery scab pathogen in suspicious samples; in particular to distinguish powdery scab from very similar symptoms of common scab caused by *Streptomyces* sp. (Fig. 2). Hence, it is a valuable tool for the disease assessment in seed certification as part of a sustainable disease management of powdery scab.

In this context, the Sss AgriStrip plays an important part in the implementation of the postulated preventive control concept: „Clean seed in clean soil“ described in Merz U. and R.E. Falloon, 2009. Review: Powdery scab of potato - Increased knowledge of pathogen biology and disease epidemiology for effective disease management. Potato Research 52 (1), 17-37.

Test principle

The rapid assay Sss AgriStrip is a lateral-flow immunochromatographic test based on an antigen-antibody reaction that is initiated by inserting the strip into the sample extract. The sample liquid migrates upwards and red colored lines will become visible within a few minutes.

Both test and control lines become visible with positive extracts (containing Sss), whereas negative samples produce only the upper control line (Fig. 1). Intense coloration is reached within 10-15 min and the result can be registered. Dried test strips can be kept as permanent records.

Specificity and sensitivity

The monoclonal antibodies used in this test were made against resting spores of *Spongospora subterranea* f. sp. *subterranea* (Sss) from potato tubers. All

isolates of Sss of several origins (France, Germany, Israel, Japan, New Zealand, Peru, USA, Scotland, and Switzerland) tested so far were detected. The sensitivity of the Sss AgriStrip is very similar to the DAS-ELISA format and allows the reliable detection of the pathogen in smallest scab lesions.


Cross-reaction has been observed with resting spores of the related *Spongospora subterranea* f. sp. *nasturtii* which infects watercress but none with other related fungi such as *Plasmodiophora brassicae*. No cross-reaction was observed with *Streptomyces* sp. (the causal agent of common scab). Hence this assay is also suitable for distinguishing lesions on potatoes caused by Sss from those caused by *Streptomyces* sp., *Rhizoctonia solani* and other fungi.

The reactivity is limited to resting spores on potato tuber or in melanized „old“ root galls. Other stages of Sss such as zoospores and zoosporangia in roots do not and developing resting spores in „young“ root galls may not react.

In contrast to the DAS-ELISA format, resting spores can not be tested in soil samples with this AgriStrip due to nonspecific interaction of soil particles with the lateral flow test.

Please read the notes on page 3 before starting the assay.

Fig. 1
Test strip

Test procedure

Sampling and extraction

Wash the tubers and cut with a sharp and disinfected knife small pieces (1-2 cm², weight approx. 0.2 g) of tuber peel (Fig. 3a) containing scab symptoms (Fig. 2a + 2b).

Place this sample in an extraction bag and add 4 ml of extraction buffer A with a disposable pipette*.

Homogenize the tissue thoroughly with the hand-held homogenizer on a flat surface (Fig. 3b).

Transfer 4 drops (approx. 150 µl) of this extract into a disposable cuvette with a disposable pipette* (Fig. 3c).

Test

Place a strip with the end marked «sample» into the extract (Fig. 3d). Observe the coloration of the lines that is completed after 15 min and can then be evaluated. Avoid any wetting of the strip above the sample mark with plant extract.

To prevent any confusion, label extraction bags, cuvettes and strips with a permanent marker.

Analysis of results

- 1) Both, a clearly visible test and control line are obtained with extracts containing high number of sporeballs as shown in Fig. 1. This means the test is **positive**.
- 2) A test line which appears weaker than the control line is obtained with extracts containing less sporeballs. This means the test is **positive**.
- 3) The test is **negative**, if the control line but no test line appears after 15 min. This means that no Sss sporeballs are in the extract or in a number below the detection limit.
- 4) If neither test nor control line become visible, the test is **invalid** and should be repeated with a fresh strip.
- 5) Rarely, a faint test line might become visible after 20 – 30 min. Repeat the test in this case with a new extract or send the sample to a laboratory for verification with another method.

Dried test strips can be kept as permanent records, even though the coloration of bands is stronger when the strips are still wet.

Fig. 2a. Symptoms of common scab




Fig. 2b. Symptoms of powdery scab

* Use of the disposable pipettes

Do not accidentally contaminate the extraction buffer. Please make sure that the disposable pipettes are immediately disinfected and disposed of, after having been in contact with samples. Use always a **fresh pipette** from the pouch for dosing extraction buffer.

Fig. 3. AgriStrip Test procedure

Notes

- 1) Store the strips and the extraction buffer at 4°C. However, limited exposure to ambient temperature (10 – 30°C), such as during transport and use in the field, does not affect the quality. Keep the packaging of the strips (containing desiccant bags) always hermetically closed. Absorbed moisture can lead to poor results or even complete failure of the test.
- 2) Strips must be used before the expiration date indicated on the label of the packaging.
- 3) Use the Sss AgriStrip always with the Extraction buffer A supplied with the strips. Other buffers might fail.
- 4) It is recommended to run the test at ambient temperature (15 – 25°C). At temperatures below 15°C, the speed of color development of test and control lines slows down.
- 5) The extraction buffer contains sodium azide as preservative. Keep out of reach of children and do not ingest.

Hygiene measures

To avoid spread of the pathogen, all disposables (extraction bags, pipettes, and cuvettes) must be properly disinfected. This can be done by submerging for at least 60 min in a freshly prepared 0.1 % bleach solution (sodium hypochlorite). After that, this material can be disposed of together with the tuber residues of the samples according to local regulations. Please contact local phytosanitary authorities for specific recommendations.

Fig. Sss AgriStrip Complete kit 25

Content of the AgriStrip Complete kit 25

- 25 strips packed with desiccant bags
- 100 ml AgriStrip Extraction buffer A
- 25 pipettes, disposable
- 25 extraction bags
- 25 cuvettes, disposable
- 1 cuvette rack

Note: the hand homogenizer is not included in the Complete kit 25.

Content of the AgriStrip Set

AgriStrip Set 25

- 25 strips packed with desiccant bags
- 100 ml AgriStrip Extraction buffer, ready-to-use A

AgriStrip Set 100

- 100 strips packed with desiccant bags
- 500 ml AgriStrip Extraction buffer, ready-to-use A

Ordering Information

Product	Art. No.	Assays
Sss AgriStrip Complete kit 25	111181	25
Sss AgriStrip Set 25	111182	25
Sss AgriStrip Set 100	111183	100
Optional Products	Art. No.	Size
Cuvette rack, holds 12 cuvettes	2166	1
Cuvettes, disposable	2534	100
Pipettes, disposable	2292	500
Extraction bags Universal	430100	100
Homogenizer hand model	400010	1

Disclaimer

This BIOREBA product is guaranteed to meet the specifications mentioned in this product information. No further warranties are given. Should this product fail for reasons other than inappropriate handling or misuse, BIOREBA AG

will replace the product free of charge or refund the purchase price.

BIOREBA AG shall not be liable for any direct or indirect, special or consequential damages of any kind resulting from the use of this product.

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