

Product Information

Consumables

ELISA plates

Product	Art. No.	Size
Nunc-Immuno plates MaxiSorp F96	520000	10
Standard packaging 6 x 10 plates	520000	60
ELISA plate covers (for multiple use)*	520001	10
Rubber sheets for ELISA plates (for multiple use)*	520003	10



* ELISA plates are covered during incubations with a sheet of Parafilm which is firmly pressed down by a rubber sheet inside a plate cover. The plate cover is fixed to the plate with rubber bands. This prevents evaporation during incubation.

Reservoirs

Product	Art. No.	Size
Single reservoir 70 ml (polypropylene)	330.01	20
8-channel reservoir 40 ml (polystyrene)	330.08.9	10
12-channel reservoir 40 ml (polystyrene)	330.12.9	10



1/1
Version: 5 - 06.12.2019
Adaptions from last version: pipettes and tips removed.