

Product Information

Consumables and Pipettes

ELISA plates

Product	Art. No.	Size
Nunc-Immuno plates MaxiSorp F96	520000	10 plates
Standard packaging: 6 x 10 plates	520000	60 plates
ELISA plate covers (for multiple use)*	520001	10 covers
Rubber sheets for ELISA plates (for multiple use)*	520003	10 sheets

* ELISA plates are covered during incubations with a sheet of Parafilm which is firmly pressed down by a rubber sheet inside a plate cover. The plate cover is fixed to the plate with rubber bands. This prevents evaporation during incubation.



Digital micro pipettes Acura® and Qualitips

Product	Art. No.	Size
Single channel pipette 0.5-10 µl	825.0010	1
Microtips (natural) 10 µl, bulk	309.0010B	1000 / bag
Microtips (natural) 10 µl	309.0010R	10 x 96 / rack
Single channel pipette 20-200 µl	825.0200	1
Pipette tips (natural) 350 µl, bulk	308.0350B	1000 / bag
Pipette tips (natural) 350 µl	308.0350R	10 x 96 / rack
Single channel pipette 100-1000 µl	825.1000	1
Pipette tips (blue) 1000 µl, bulk	319.1000B	4 x 250 / bag
Pipette tips (blue) 1000 µl	319.1000R	8 x 60 / rack



Digital multichannel pipette Acura®

Product	Art. No.	Size
8-channel pipette 40-350 µl	855.08.350	1
12-channel pipette 40-350 µl	855.12.350	1



Reservoirs

Product	Art. No.	Size
Single reservoir 70 ml (polypropylene)	330.01	20
8-channel reservoir 40 ml (polystyrene)	330.08.9	10
12-channel reservoir 40 ml (polystyrene)	330.12.9	10

