

BIOrack

The perfect partner for BIOREBA Extraction Bags

A CALIBRE SCIENTIFIC COMPANY

BIOREBA



- **Tailored Fit:** Specifically designed for use with BIOREBA Universal extraction bags.
- **Contamination Prevention:** Fully contains any potential liquid leakage.
- **High Capacity:** Holds 10 or more bags simultaneously.
- **Space-Saving:** Stackable design for compact storage in limited lab spaces.
- **Easy to Clean:** Washable at temperatures up to 100°C (not autoclavable).
- **Durable & Lab-Ready:** Made from ASA (Acrylonitrile Styrene Acrylate), which withstands exposure to many common lab chemicals*, offers excellent stability under UV light and remains impact-resistant at room temperature.
- **Lightweight & Portable:** Easy to move.
- **Available options:** Supplied in packs of five, in different colors

*Not resistant to acetone or chloroform

Your Partner in Agro-Diagnostics