

# Product Information

## Consumables

### ELISA plates

Product	Art. No.	Size
Nunc-Immuno plates MaxiSorp F96	520000	10
Standard packaging 6 x 10 plates	520000	60
ELISA plate covers (for multiple use)*	520001	10
Rubber sheets for ELISA plates (for multiple use)*	500003	10



\* ELISA plates are covered during incubations with a sheet of Parafilm which is firmly pressed down by a rubber sheet inside a plate cover. The plate cover is fixed to the plate with rubber bands. This prevents evaporation during incubation.

### Reservoirs

Product	Art. No.	Size
Single reservoir 70 ml (polypropylene)	330.01	20
8-channel reservoir 40 ml (polystyrene)	330.08.9	10
12-channel reservoir 40 ml (polystyrene)	330.12.9	10



1/1  
Version: 6 - 21.05.2021  
Adaptions from last version: catalog number of rubber sheets.