



▶ HOMEX 7 Video

## HOMEX 7

The HOMEX 7 is the newest homogenization device member of the BIOREBA HOMEX family. We improved it in every way: more compact and lighter compared to the previous model. It's also more efficient at homogenizing up to 5 samples using BIOREBA "Universal, Universal-U-form" extraction bags and up to 4 samples using BIOREBA "Universal-Long" extraction bags.

- ✓ Ergonomic handling
- ✓ High throughput (ca. 250 samples/h)
- ✓ Optimized for BIOREBA extraction bags
- ✓ Improved materials against corrosion
- ✓ 27 % lighter than the previous model
- ✓ 100 % Made in Switzerland

Product name	HOMEX 7
Application	Plant samples homogenizer
Power	max. 160 W
Power supply unit	100 – 240 V AC, 50/60 Hz, 2A
Max. Speed	3000 rpm / adjustable (scale 1-10)
Sound level	< 72.5 dB
Dimensions	834 mm x 451 mm x 235 mm
Weight	24.8 kg

### Ordering information

400007S	HOMEX 7 incl. Standard Rack
4000074	HOMEX 7 incl. Rack 4/25
400014S	Standard Rack for HOMEX 7
4000164	Rack 4/25 for HOMEX 7