

How to store products?



A CALIBRE SCIENTIFIC COMPANY

BIOREBA

- ELISA Coating Buffer Tablets
- Microtiter Plates
- ELISA Washing Buffer Powder and Tablets
- Potato DNA/RNA-Rapid Extraction Set

+20 °C
(Room Temperature)

- AgriStrip Kits and Sets
- ELISA Conjugate Buffer
- Egg Albumin
- ELISA Extraction Buffers
- IgG/Conjugate Antibodies
- ELISA Negative Control
- ELISA Positive Control
- ELISA Substrate Buffer

+4 °C

- qPCR Set Components
- ELISA Substrate (pNPP)

-20 °C

Your Partner in Agro-Diagnostics