



Products and Services


2023



PRODUCTS



EXTRACTION BAGS



HOMEX



POTATO DNA/RNA RAPID EXTRACTION



ELISA REAGENTS, SETS & KITS



SWISS ACCREDITATION
TESTING SERVICE



qPCR



PLANT PATHOGENS



CROPS

SWISS TESTING LABS
Association of Swiss Laboratories
Verband Schweizer Laboratorien
Association des Laboratoires Suisses
Associazione dei Laboratori Svizzeri



AGRISTRIP

The Company



BIOREBA AG was the **first** company in the world to produce and commercialise ELISA reagents for plant virus diagnostics in 1980. It is an independent company that was founded in close collaboration with Agroscope, the Swiss Centre of Competence for Agricultural and Food Research. Since then, BIOREBA has established collaborations with numerous other research institutes and universities in many countries. Since its foundation almost four decades ago, the company has grown steadily and established a reputation for **high quality diagnostics**. Today, the research community and diagnostic laboratories around the world rely on our products to potatoes, vines, fruit trees, small fruits, vegetables, ornamentals and field crops.

Our philosophy is to provide our customers with high quality products and services for their specific diagnostic needs. Customers **benefit** from our expertise and many years of experience in the development and introduction of complete assays for the detection of plant pathogens. BIOREBA has well-equipped facilities for R&D activities and production of these test systems, as well as an **accredited** testing laboratory to provide customers with serological and molecular analyses of plant pathogens. The company also provides specialised equipment and disposables for plant sample preparation and liquid handling.

BIOREBA products are available **throughout the world**. The sales and marketing department and an extensive network of distributors ensure a rapid response to customer enquiries and prompt delivery of products. Our well trained staff are committed to maintaining high quality standards.

BIOREBA has been **awarded** the internationally recognized **ISO 9011** Quality Management System certificate. In addition, our testing services laboratory is accredited according to **ISO/IEC 17025** for the analysis of plant pathogens. Both certificates are an incentive for us to continue serving you with the best quality products and to help you grow healthy crops and yields.



Please feel free to contact us. We will be happy to help you. We will usually get back to you within 24 hours:

BIOREBA AG

Christoph Merian-Ring 7
4153 Reinach
Switzerland

Phone: +41 61 712 11 25
Fax: +41 61 712 11 17
Email: admin@bioreba.ch
Web: www.bioreba.com



Sample Preparation


Page 6

Extraction Bags



Page 7 – 8

Extraction Devices



Page 9 – 10

Extraction Buffers



Plant Pathogen Diagnostic Tests

Page 14 – 16

AgriStrip – *Lateral flow test*



Page 17 – 43

ELISA



Page 44 – 46

qPCR



Additional Products

Page 47

Buffers, Substrates



Page 48

Chemicals, Consumables



Services

Page 49

Testing Service



Page 50

Quality Management



Page 51

Distribution, Mission



Extraction Bags



Our homogenization system includes **four different extraction bags** in different sizes, with or without a synthetic intermediate layer. All extraction bags are ideal for tissue extraction for use in ELISA and PCR.


The characteristics of our extraction bags are:

- Optimal sampling and sample preparation (samples can be labeled and stored)
- Optimal homogenization
- Optimal filtration of the plant extracts (bags with synthetic intermediate layer)
- Optimal recovery of the plant sap
- No cross-contamination between test samples
- Optimal for homogenizing with our HOMEX 7, HOMEX XS or Homogenizer hand model
- Suitable for homogenizing leaves, sprouts, roots, seeds, woody samples, sieved soil

Ordering information for BIORBEA Extraction bags: Each unit contains **100 bags**. For 1000 bags (one full box) you need to order the extraction bag item 10 times. i

Art. No.	Product	Quantity
430100	Extraction bags Universal 12 x 15 cm for our standard rack <ul style="list-style-type: none"> • with synthetic intermediate layer for optimal filtration, heavy duty plastic for testing all kind of samples • max. filling volume for use with HOMEX 7: 10 ml 	100
440100	Extraction bags Universal-U-form 12 x 15 cm for our standard rack <ul style="list-style-type: none"> • heavy duty plastic with U-formed compartment with synthetic intermediate layer for testing seeds • max. filling volume for use with HOMEX 7: 4 ml 	100
470100	Extraction bags Universal-Long 15 x 28 cm for our rack 4/25 <ul style="list-style-type: none"> • heavy duty plastic with 15 mm wide sealing borders on bottom and sides with synthetic intermediate layer • max. filling volume for use with HOMEX 7: 50 ml 	100
420100	Extraction bags Long 15 x 40 cm <ul style="list-style-type: none"> • heavy duty plastic without intermediate layer for testing long leaves • max. filling volume for use with HOMEX 7: 200 ml 	100

Extraction Devices



Our **homogenizers** are ideal for tissue extraction of samples for use in ELISA and PCR. Samples remain in extraction bags throughout sample processing, eliminating cross-contamination and the need for washing between samples.

HOMEX 7 incl. Standard Rack – *The best-selling homogenizer*


- Semi-automated homogenizer
- Ergonomic and easy handling
- Tissue samples of different thickness may be homogenized
- High throughput (ca. 250 samples/h can be processed)
- Can be equipped with two different racks (see Rack 4/25) to hold different sizes of BIOREBA extraction bags
- Changing the rack is fast and easy
- 27% lighter than the previous model
- 100% Made in Switzerland

HOMEX 7 incl. Rack 4/25

- Semi-automated HOMEX 7 with a optional longer rack
- Holds 1-4 BIOREBA extraction bags Universal-Long
- High throughput (ca. 250 samples/h can be processed)
- Tissue samples of different thickness may be homogenized
- Changing the rack is fast and easy
- Ergonomic and easy handling
- 27% lighter than the previous model
- 100% Made in Switzerland



Scan me.
Watch the HOMEX video
on our YouTube channel.



HOMEX XS – *The compact homogenizer*

- Automated easy-to-use system
- Space saving and low weight machine
- High efficiency with homogenization of more than 100 samples per hour
- For homogenization of all kinds of tissues without sharp edges
- Filtered extract can be used for analyses such as ELISA, PCR or other analytical applications
- Optimized for BIOREBA extraction bags Universal and Universal-U-form
- Adjustable rotation speed and homogenization time
- Stop timer function
- Ergonomic handling

Extraction Devices



Homogenizer hand model

- Suitable for small sample numbers and for samples of softer tissue (leaves, sprouts, flowers etc.)
- Optimal grinding
- Can be used with all types of BIORÉBA extraction bags

Products	Art. No.	Quantity
HOMEX 7 incl. Standard Rack	400007S	1
HOMEX 7 incl. Rack 4/25	4000074	1
HOMEX XS	400008	1
Homogenizer hand model	400010	1

Replacement parts for HOMEX 7	Art. No.	Quantity
Standard Rack for HOMEX 7	400014S	1
Rack 4/25 for HOMEX 7	4000164	1
Replacement head for HOMEX 7	400017	1
Aluminium bars for Standard Rack HOMEX 7	400027	6



BIOpick is a device for taking heel core samples from potato tubers (e.g. for testing for quarantine bacterial diseases such as ring rot caused by *Clavibacter* and brown rot caused by *Ralstonia*). It makes coring more ergonomic and allows up to 200 heel end cores to be collected and transferred into a 50ml collection tube. Sampling can be performed manually or semi-automatically by connecting the BIOpick to a laboratory blender using the supplied adapter.

BIOpick – Your Potato Sample Device

- Time-saving
- Ergonomic
- Reusable
- Flexible use: Manual or Semi-automatic
- High throughput: up to 20 cores/min



Scan me.
Watch the BIOpick video on our YouTube channel.



BIOpick – Your Potato Sample Device	Art. No.	Quantity
BIOpick, stainless steel, 3pcs	BP003	1
Adapter for BIOpick, stainless steel	BPA01	1
Brushes for BIOpick, 3 pcs	BPB03	1 set
BIOpick set (BIOpick, adapter, cleaning brushes)	BPS01	1 set

Potato DNA/RNA rapid extraction set



Nucleic acid extraction for potato samples

The Potato DNA/RNA Rapid Extraction Kit is designed for phytodiagnostic use. It allows **quick and easy** extraction of DNA and RNA from potato tuber and leaf samples. The extraction requires the use of only two buffers. Time-consuming processing steps are eliminated. A particularly cost-effective advantage is that the Potato DNA / RNA Rapid Extraction Kit does not require expensive machinery.

The characteristics of our Potato DNA/RNA rapid extraction set are:


- Fast extraction in 30 minutes only
- Easy to use: only 5 steps
- Optimized for pools of up to 25 potato samples
- No specialised machines are needed
- Room temperature storage and shipment
- Total nucleic acids – extracts DNA and RNA in one procedure
- Validated and accredited method according to ISO/IEC 17025
- Extracted nucleic acids can be used for diagnostic tests such as real time PCR and PCR macroarray
- Assays for viruses, bacteria and fungi

Sets are stored and shipped at **room temperature** for easy handling.

The set contains all the reagents and extraction bags required for DNA/RNA extraction. Separate samples for DNA and RNA extraction are not required as the kit extracts total nucleic acids. Downstream applications such as real-time PCR or PCR macroarray for detection of viruses, bacteria or fungi can be performed using one extract.

The extraction method using the Potato DNA/RNA Rapid Extraction Kit has been validated and accredited according to ISO / IEC 17025. The method is **approved for seed potato certification**.

Four set sizes are available for pools of 10 or up to 25 potato samples and for 96 or 192 extractions respectively (see next page).



Art. No.	Product	No. of reactions	Amount of samples (potatoes)
810003	Potato DNA/RNA rapid extraction set 96 / 10x for 96 extractions (Pool size: up to 10 tuber samples) set includes: <ul style="list-style-type: none"> • Rapid extraction buffer EB 1 (p) 100ml • Rapid extraction buffer EB 2 (p) 50 ml • Extraction bags Standard 12x15 cm 100 pcs 	96	up to 10
810004	Potato DNA/RNA rapid extraction set 96 / 25x for 96 extractions (Pool size: up to 25 tuber samples) set includes: <ul style="list-style-type: none"> • Rapid extraction buffer EB 1 (p) 200ml • Rapid extraction buffer EB 2 (p) 50 ml • Extraction bags Standard 12x15 cm 100 pcs 	96	up to 25
810005	Potato DNA/RNA rapid extraction set 192 / 10x for 192 extractions (Pool size: up to 10 tuber samples) set includes: <ul style="list-style-type: none"> • Rapid extraction buffer EB 1 (p) 200ml • Rapid extraction buffer EB 2 (p) 100 ml • Extraction bags Standard 12x15 cm 2 x 100 pcs 	192	up to 10
810006	Potato DNA/RNA rapid extraction set 192 / 25x for 192 extractions (Pool size: up to 25 tuber samples) set includes: <ul style="list-style-type: none"> • Rapid extraction buffer EB 1 (p) 2 x 200ml • Rapid extraction buffer EB 2 (p) 100 ml • Extraction bags Standard 12x15 cm 2 x 100 pcs 	192	up to 25

Extraction buffers for ELISA



ELISA Extraction buffer

Depending on the plant samples to be analyzed, a **different** extraction buffer will be required. For most plants, the extraction buffer «General» is suitable. For the analysis of grapevines, the extraction buffer «Grapevine» is required. For testing tubers and sprouts of potatoes, bulbs of onions and other *Allium* spp. and seeds, the extraction buffer «Bulbs and Tubers» is required. For CaMV tests the special extraction buffer «Caulimo» is required. Information on the correct use of the extraction buffer is given in the «Product Information» supplied with our ELISA reagents.

Product	Art. No.	Size	sufficient for
Extraction buffer «General» (10x)	110120	500 ml	1000 samples**
Extraction buffer «Grapevine» (5x)	110123	2 x 500 ml	1000 samples**
Extraction buffer «Bulbs & Tubers» ***	110122	500 ml	1000 samples**
Extraction buffer «Caulimo» ****	110124	500 ml	1000 samples
Extraction buffer «Caulimo» ****	110224	250 ml	500 samples

** Calculated with 5 ml extraction buffer per sample of 0.25 g

*** Composed of 500 ml extraction buffer «General» (10x) and 50 g egg albumin

**** Includes extraction buffer «Caulimo» (10x), urea, and egg albumin

Overview: Plant Pathogen Diagnostic Tests


Our **AgriStrip** products are listed in alphabetical order on page 15. Our **ELISA** reagents, reagent sets and complete kits are listed in alphabetical order on page 17–42. Our **qPCR** products are listed on page 43–45. **The table below illustrates the prevalence of pathogens in relation to different crop groups:**

Pathogen name	Acronym	AgriStrip	ELISA	qPCR
Potatoes				
Alfalfa mosaic virus	AMV		•	
Andean potato latent virus	APLV		•	
Andean potato mottle virus	APMoV		•	
Potato leafroll virus	PLRV		•	•
Potato mop-top virus	PMTV		•	
Poty group	Poty group		•	
Potato virus A	PVA	•	•	•
Potato virus M	PVM	•	•	•
Potato virus S	PVS	•	•	
Potato virus V	PVV		•	
Potato virus X	PVX	•	•	
Potato virus Y (monoclonal cocktail)	PVY-mono cock	•	•	•
Potato virus Y (monoclonal)	PVY-mono	•	•	•
Potato virus Y (polyclonal)	PVY-poly	•	•	•
Potato virus Y-necrotic	PVY-necrotic	•	•	•
<i>Ralstonia solanacearum</i>	Rs	•	•	
<i>Spongospora subterranea f.sp. subterranea</i>	Sss	•	•	
Tomato black ring virus	TBRV		•	
Tobacco necrosis virus	TNV		•	
Tomato spotted wilt virus	TSWV	•	•	
<i>Verticillium spp.</i>	Verticillium		•	
Grapevines				
Arabis mosaic virus	ArMV		•	
ArMV+GFLV	ArMV+GFLV		•	
Grapevine fanleaf virus	GFLV		•	
Grapevine fleck virus	GFkV		•	
Grapevine leafroll assoc. virus 1	GLRaV-1		•	
Grapevine leafroll assoc. virus 1+3	GLRaV-1+3		•	
Grapevine leafroll assoc. virus 2	GLRaV-2		•	
Grapevine leafroll assoc. virus 3	GLRaV-3		•	
Grapevine leafroll assoc. virus generic 4 strains	GLRaV-4 strains		•	
Grapevine leafroll assoc. virus 4 strain 6	GLRaV-4 strain 6		•	
Grapevine pinot gris virus	GPGV		•	
Grapevine virus A	GVA		•	
Grapevine red blotch virus	GRBV			•
Raspberry ringspot virus-ch	RpRSV-ch		•	
Raspberry ringspot virus-g	RpRSV-g		•	
Strawberry latent ringspot virus	SLRSV		•	
Tomato black ring virus	TBRV		•	
Tomato ringspot virus	ToRSV		•	
Tomato ringspot virus-Ch	ToRSV-Ch		•	
Tobacco ringspot virus	TRSV		•	
Vegetables				
<i>Acidovorax avenae subsp. citrulli</i>	Aac	•	•	
African cassava mosaic virus	ACMV (TYLCV)		•	
Bean common mosaic virus	BCMV		•	
Bean common mosaic necrosis virus	BCMNV		•	
Bean golden mosaic virus	BGMV (TYLCV)		•	
Beet necrotic yellow vein virus	BNYVV		•	
Cauliflower mosaic virus	CaMV		•	
Cucumber green mottle mosaic virus	CGMMV		•	
Cucumber mosaic virus	CMV	•	•	

Pathogen name	Acronym	AgriStrip	ELISA	qPCR
Vegetables				
Cucurbit yellow stunting disorder virus	CYSDV		•	
Erysimum latent virus	ErLV		•	
Garlic common latent virus	GCLV		•	
Lettuce mosaic virus	LMV	•	•	
Leek yellow stripe virus	LYSV		•	
Mirafiori lettuce big-vein virus	MiLBVV		•	
Onion yellow dwarf virus	OYDV		•	
Pepino mosaic virus	PepMV	•	•	
Pepper mild mottle virus	PMMoV	•	•	
Poty group	Poty group		•	
Papaya ringspot virus	PRSV (=WMV-1)		•	
Potato virus X	PVX	•	•	
Potato virus Y	PVY	•	•	
Radish mosaic virus	RaMV		•	
Soybean mosaic virus	SMV		•	
Squash leaf curl virus	SLCV		•	
Squash mosaic virus	SqMV		•	
Tomato bushy stunt virus	TBSV		•	
Tomato chlorotic spot virus	TCSV (TSWV)		•	
Tobacco mosaic virus	TMV	•	•	
Tobacco necrosis virus	TNV		•	
Tomato apex necrosis virus	ToANV		•	
Tomato brown rugose fruit virus	ToBRFV	•	•	•
Tomato mosaic virus	ToMV	•	•	
Tomato mottle virus	ToMoV (TYLCV)		•	
Tomato ringspot virus	ToRSV		•	
Tomato ringspot virus-Ch	ToRSV-Ch		•	
Tospovirus broad-spectrum (I,II,III)	Tospo (I,II,III)		•	
Tobacco streak virus	TSV		•	
Tomato spotted wilt virus	TSWV	•	•	
Turnip mosaic virus	TuMV		•	
Tomato yellow leaf curl virus	TYLCV		•	
Turnip yellow mosaic virus	TYMV		•	
<i>Verticillium spp.</i>	Verticillium		•	
Watermelon mosaic virus 1	WMV-1 (=PRSV)		•	
Watermelon mosaic virus 2	WMV-2		•	
Zucchini yellow mosaic virus	ZYMV		•	
Field Crops				
Alfalfa mosaic virus	AMV		•	
Barley yellow dwarf virus-B	BYDV-B		•	
Barley yellow dwarf virus-F	BYDV-F		•	
Barley yellow dwarf virus-RPV	BYDV-RPV		•	
Beet necrotic yellow vein virus	BNYVV		•	
Cauliflower mosaic virus	CaMV		•	
Cucumber mosaic virus	CMV	•	•	
Cereal yellow dwarf virus-RPV	CYDV-RPV (BYDV-RPV)		•	
Groundnut ringspot virus	GRSV (TSWV)		•	
Maize chlorotic mottle virus	MCMV	•	•	
Maize dwarf mosaic virus	MDMV		•	
Petunia asteroid mosaic virus	PeAMV		•	
Poty group	Poty group		•	
Potato virus X	PVX	•	•	
Potato virus Y	PVY	•	•	
Radish mosaic virus	RaMV		•	
Sugarcane mosaic virus	SCMV		•	
Tobacco mosaic virus	TMV	•	•	
Tobacco necrosis virus	TNV		•	
Tobacco streak virus	TSV		•	
Tomato spotted wilt virus	TSWV	•	•	
Turnip yellow mosaic virus	TYMV		•	

Pathogen name	Acronym	AgriStrip	ELISA	qPCR
Ornamentals				
<i>Acidovorax avenae subsp. cattleyae</i>	Aacat		•	
Alfalfa mosaic virus	AMV		•	
Apple mosaic virus	ApMV		•	
Arabis mosaic virus	ArMV		•	
Calibrachoa mottle virus	CbMV		•	
Cucumber mosaic virus	CMV	•	•	
Erysimum latent virus	ErLV		•	
Impatiens necrotic spot virus	INSV	•	•	
Petunia asteroid mosaic virus	PeAMV		•	
Pelargonium flower-break virus	PFBV		•	
Pelargonium leaf curl virus	PLCV		•	
Pelargonium line pattern virus	PLPV		•	
Prunus necrotic ringspot virus	PNRSV		•	
Poty group	Poty group		•	
Potato virus A	PVA	•	•	
Potato virus X	PVX	•	•	
Potato virus Y	PVY	•	•	
<i>Ralstonia solanacearum</i>	Rs	•	•	
Strawberry latent ringspot virus	SLRSV		•	
Tobacco mosaic virus	TMV	•	•	
Tobacco necrosis virus	TNV		•	
Tomato mosaic virus	ToMV		•	
Tomato ringspot virus	ToRSV		•	
Tomato ringspot virus-Ch	ToRSV-Ch		•	
Tospovirus broad-spectrum (I,II,III)	Tospo (I,II,III)		•	
Tobacco ringspot virus	TRSV		•	
Tobacco streak virus	TSV		•	
Tomato spotted wilt virus	TSWV	•	•	
<i>Verticillium spp.</i>	Verticillium		•	
<i>Xanthom. campestris pv. pelargonii</i>	Xcp (Xh)		•	
Fruit Trees & Small Fruits				
Apple chlorotic leaf spot virus	ACLSV		•	
Apple mosaic virus	ApMV		•	
Apple proliferation phytoplasma	ApP		•	
Arabis mosaic virus	ArMV		•	
Apple stem grooving virus	ASGV		•	
Apple stem pitting virus	ASPV		•	
Cherry leaf roll virus	CLRV		•	
Citrus tristeza virus	CTV		•	
<i>Erwinia amylovora</i>	Ea	•		
Petunia asteroid mosaic virus	PeAMV		•	
Papaya ringspot virus	PRSV (=WMV-1)		•	
Plum pox virus (Sharka)	PPV	•	•	
Prune dwarf virus	PDV		•	
Prunus necrotic ringspot virus	PNRSV		•	
Poty group	Poty group		•	
Raspberry bushy dwarf virus	RBDV		•	
Raspberry ringspot virus-ch	RpRSV-ch		•	
Soybean mosaic virus	SMV		•	
Strawberry latent ringspot virus	SLRSV		•	
Strawberry mild yellow edge PV	SMYEPV		•	
Tomato ringspot virus-ch	ToRSV-ch		•	
Tomato ringspot virus	ToRSV		•	
Tobacco ringspot virus	TRSV		•	
Tobacco streak virus	TSV		•	
<i>Verticillium spp.</i>	Verticillium		•	
<i>Xanthomonas fragariae</i>	Xfr		•	

AgriStrip – Lateral flow test



The rapid one-step assay for the detection of plant pathogens

The **AgriStrip**, developed and manufactured by BIOREBA, is based on lateral flow immunochromatography using specific antibodies. The AgriStrip is designed to confirm the presence of a specific pathogen in suspect plant samples on the spot **within minutes**.

The characteristics of our AgriStrip are:

- The performance of the test is very easy and fast
- Cost effective
- No need of special laboratory equipment
- Especially suitable when testing is required on-site, such as in greenhouses or fields
- Test strips can be kept as permanent records
- Available as Complete kit for 25 assays or as set with 25 or 100 strips, respectively

AgriStrip Complete Kit 25

Component	Assays	
	25	
AgriStrips	25 pcs	
AgriStrip Extraction buffer ready-to-use	100 ml	
Pipettes, disposable	25 pcs	
Extraction bags Universal	25 pcs	
Cuvettes, disposable	25 pcs	
Cuvette rack	1 pc	

AgriStrip Set

Component	Assays	
	25	100
AgriStrips	25 pcs	100 pcs
AgriStrip Extraction buffer ready-to-use	100 ml	500 ml

Optional Products

Art. No.		Size
2166	Cuvette rack, holds 12 cuvettes	1
2534	Cuvettes, disposable	100
2292	Pipettes, disposable	500
430100	Extraction bags Universal	100
400010	Homogenizer hand model	1

MOQ: to further increase our production efficiency, we have decided to offer some of our AgriStrips (Lateral Flow Tests) on a minimum order quantity (MOQ) basis. The following pathogens will be affected: Aac, INSV, PMMoV, PVA, PVM, PVS, PVX and Sss. The MOQ only applies to the 8 plant pathogens mentioned, other pathogens such as CMV, EA, or PPV for example are not affected and have (as before) no MOQ.



Pathogen	Acronym	Art. No.	Format	Assays
<i>Acidovorax avenae</i> subsp. <i>citrulli</i> (on request)	Aac	161981	AgriStrip Complete kit 25	25
		161982	AgriStrip Set 25	25
		161983	AgriStrip Set 100	100
Cucumber mosaic virus	CMV	160681	AgriStrip Complete kit 25	25
		160682	AgriStrip Set 25	25
		160683	AgriStrip Set 100	100
<i>Erwinia amylovora</i> (Fire blight)	Ea	153081	AgriStrip Complete kit 25	25
		153082	AgriStrip Set 25	25
		153083	AgriStrip Set 100	100
Impatiens necrotic spot virus (on request)	INSV	190281	AgriStrip Complete kit 25	25
		190282	AgriStrip Set 25	25
		190283	AgriStrip Set 100	100
Lettuce mosaic virus	LMV	160581	AgriStrip Complete kit 25	25
		160582	AgriStrip Set 25	25
		160583	AgriStrip Set 100	100
Maize chlorotic mottle virus	MCMV	140781	AgriStrip Complete kit 25	25
		140782	AgriStrip Set 25	25
		140783	AgriStrip Set 100	100
Pepino mosaic virus	PepMV	161581	AgriStrip Complete kit 25	25
		161582	AgriStrip Set 25	25
		161583	AgriStrip Set 100	100
Pepper mild mottle virus (on request)	PMMoV	161781	AgriStrip Complete kit 25	25
		161782	AgriStrip Set 25	25
		161783	AgriStrip Set 100	100
Plum pox virus (Sharka)	PPV	150581	AgriStrip Complete kit 25	25
		150582	AgriStrip Set 25	25
		150583	AgriStrip Set 100	100
Potato virus A (on request)	PVA	112181	AgriStrip Complete kit 25	25
		112182	AgriStrip Set 25	25
		112183	AgriStrip Set 100	100
Potato virus M (on request)	PVM	112281	AgriStrip Complete kit 25	25
		112282	AgriStrip Set 25	25
		112283	AgriStrip Set 100	100
Potato virus S	PVS	110381	AgriStrip Complete kit 25	25
		110382	AgriStrip Set 25	25
		110383	AgriStrip Set 100	100
Potato virus X (on request)	PVX	110481	AgriStrip Complete kit 25	25
		110482	AgriStrip Set 25	25
		110483	AgriStrip Set 100	100
Potato virus Y	PVY	112981	AgriStrip Complete kit 25	25
		112982	AgriStrip Set 25	25
		112983	AgriStrip Set 100	100
<i>Ralstonia solanacearum</i>	Rs	170581	AgriStrip Complete kit 25	25
		170582	AgriStrip Set 25	25
		170583	AgriStrip Set 100	100
<i>Spongospora subterranea</i> f.sp. <i>subterranea</i> (Powdery scab of potato) (on request)	Sss	111181	AgriStrip Complete kit 25	25
		111182	AgriStrip Set 25	25
		111183	AgriStrip Set 100	100
Tobacco mosaic virus	TMV	190481	AgriStrip Complete kit 25	25
		190482	AgriStrip Set 25	25
		190483	AgriStrip Set 100	100
Tomato mosaic virus	ToMV	152681	AgriStrip Complete kit 25	25
		152682	AgriStrip Set 25	25
		152683	AgriStrip Set 100	100
Tomato spotted wilt virus	TSWV	190181	AgriStrip Complete kit 25	25
		190182	AgriStrip Set 25	25
		190183	AgriStrip Set 100	100

AgriStrip^{cassette} – The rapid one-step assay



The rapid one-step assay for the detection of plant pathogens in a new format

The «**ToBRFV AgriStrip^{cassette}**» rapid test, developed and produced by BIOREBA, allows the presence of ToBRFV, also known as «Jordan virus», to be detected in suspect samples in the field **within minutes**.


ToBRFV belongs to the genus Tobamovirus and causes various symptoms such as yellow or brown wrinkled spots on fruit, narrowed and mosaic discoloration of leaves, and deformed or irregularly ripened fruit. The main host plants of ToBRFV are tomato (*Solanum lycopersicum*) and pepper (*Capsicum spp.*). ToBRFV is on the EPPO-A2 list of pests recommended for regulation as a quarantine pest as of September 2020.

The content of AgriStrip^{cassette} kit:

- 5 cassettes packed with desiccant bags
- 5 tubes containing glass beads and extraction buffer
- 5 dropper tips
- 1 technical information
- 1 product information

Pathogen	Acronym	Art. No.	Format	Assays
Tomato brown rugose fruit virus	ToBRFV	152881NF	ToBRFV AgriStrip ^{cassette} kit	5

ELISA Reagents, Sets and Kits



Individual ELISA Reagents

IgG and/or conjugate for the working volume of 200 µl/test/well available for 96, 480 or 960 assays.



ELISA Reagent Sets

Component	Assays		
	480	960	5000 [®]
IgG	0.1 ml	0.2 ml	1 ml
Conjugate	0.1 ml	0.2 ml	1 ml
Positive control	2.5 ml	2.5 ml	5 x 2.5 ml
Negative control	2.5 ml	2.5 ml	5 x 2.5 ml
Microtiter plates	5	10	60

[®] Available for potato viruses PLRV, PVA, PVM, PVS, PVX and PVY.



ELISA Complete Kits

Component	Assays			
	96	480	960	5000 [®]
IgG	0.025 ml	0.1 ml	0.2 ml	1 ml
Conjugate	0.025 ml	0.1 ml	0.2 ml	1 ml
Positive control	2.5 ml	2.5 ml	2.5 ml	5 x 2.5 ml
Negative control	2.5 ml	2.5 ml	2.5 ml	5 x 2.5 ml
Extraction buffer (10x)*	50 ml	250 ml	500 ml	5 x 500 ml
Coating buffer	1 tablet	1 tablet	2 tablets	10 tablets
Conjugate buffer (10x)	2 ml	10 ml	20 ml	100 ml
Substrate buffer (5x)	4 ml	20 ml	2 x 20 ml	200 ml
Washing buffer	1 tablet	1 pouch	2 pouches	2 containers
Substrate (pNPP)**	4 tablets (5 mg)	5 tablets (20 mg)	10 tablets (20 mg)	bulk powder (5 g)
Egg Albumin**	5 g	25 g	50 g	250 g
Microtiter plates	1	5	10	60

* For grapevine testing, twice as much buffer is included since the concentration is 5x

** Is included only with the extraction buffer «Bulbs & Tubers»

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Acidovorax avenae subsp. citrulli polyclonal - polyclonal*	Aac	161912	IgG	0.2 ml	1000
		161915	IgG	0.1 ml	500
		161919	IgG	0.025 ml	100
		161922	Conjugate	0.2 ml	1000
		161925	Conjugate	0.1 ml	500
		161929	Conjugate	0.025 ml	100
		161953	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161962	Reagent set 960		960
		161965	Reagent set 480		480
		161972	Complete kit 960		960
		161975	Complete kit 480		480
		161977	Complete kit 96		96
		Acidovorax avenae subsp. cattleyae polyclonal - polyclonal	Aacat	170612	IgG
170615	IgG			0.1 ml	500
170619	IgG			0.025 ml	100
170622	Conjugate			0.2 ml	1000
170625	Conjugate			0.1 ml	500
170629	Conjugate			0.025 ml	100
170653	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170662	Reagent set 960				960
170665	Reagent set 480				480
170672	Complete kit 960				960
170675	Complete kit 480				480
170677	Complete kit 96				96
African cassava mosaic virus	ACMV			see Tomato yellow leaf curl virus	
Alfalfa mosaic virus polyclonal - polyclonal	AMV	140512	IgG	0.2 ml	1000
		140515	IgG	0.1 ml	500
		140519	IgG	0.025 ml	100
		140522	Conjugate	0.2 ml	1000
		140525	Conjugate	0.1 ml	500
		140529	Conjugate	0.025 ml	100
		140553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		140562	Reagent set 960		960
		140565	Reagent set 480		480
		140572	Complete kit 960		960
		140575	Complete kit 480		480
		140577	Complete kit 96		96
		Andean potato latent virus polyclonal - polyclonal	APLV	110912	IgG
110915	IgG			0.1 ml	500
110919	IgG			0.025 ml	100
110922	Conjugate			0.2 ml	1000
110925	Conjugate			0.1 ml	500
110929	Conjugate			0.025 ml	100
110953	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
110962	Reagent set 960				960
110965	Reagent set 480				480
110972	Complete kit 960				960
110975	Complete kit 480				480
110977	Complete kit 96				96

*refers to coating IgG and conjugate, respectively.
Product Information on www.bioreba.com

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Andean potato mottle virus polyclonal - polyclonal	APMoV	111012	IgG	0.2 ml	1000
		111015	IgG	0.1 ml	500
		111019	IgG	0.025 ml	100
		111022	Conjugate	0.2 ml	1000
		111025	Conjugate	0.1 ml	500
		111029	Conjugate	0.025 ml	100
		111053	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		111062	Reagent set 960		960
		111065	Reagent set 480		480
		111072	Complete kit 960		960
		111075	Complete kit 480		480
		111077	Complete kit 96		96
		Apple chlorotic leaf spot virus polyclonal - polyclonal	ACLSV	151012	IgG
151015	IgG			0.1 ml	500
151019	IgG			0.025 ml	100
151022	Conjugate			0.2 ml	1000
151025	Conjugate			0.1 ml	500
151029	Conjugate			0.025 ml	100
151053	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
151062	Reagent set 960				960
151065	Reagent set 480				480
151072	Complete kit 960				960
151075	Complete kit 480				480
151077	Complete kit 96				96
Apple mosaic virus polyclonal - polyclonal	ApMV			150712	IgG
		150715	IgG	0.1 ml	500
		150719	IgG	0.025 ml	100
		150722	Conjugate	0.2 ml	1000
		150725	Conjugate	0.1 ml	500
		150729	Conjugate	0.025 ml	100
		150753	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		150762	Reagent set 960		960
		150765	Reagent set 480		480
		150772	Complete kit 960		960
		150775	Complete kit 480		480
		150777	Complete kit 96		96
		Apple proliferation phytoplasma monoclonal - monoclonal	ApP	151712	IgG
151715	IgG			0.1 ml	500
151719	IgG			0.025 ml	100
151722	Conjugate			0.2 ml	1000
151725	Conjugate			0.1 ml	500
151729	Conjugate			0.025 ml	100
151753	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
151762	Reagent set 960				960
151765	Reagent set 480				480
151772	Complete kit 960				960
151775	Complete kit 480				480
151777	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Apple stem grooving virus polyclonal - polyclonal	ASGV	150812	IgG	0.2 ml	1000
		150815	IgG	0.1 ml	500
		150819	IgG	0.025 ml	100
		150822	Conjugate	0.2 ml	1000
		150825	Conjugate	0.1 ml	500
		150829	Conjugate	0.025 ml	100
		150853	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		150862	Reagent set 960		960
		150865	Reagent set 480		480
		150872	Complete kit 960		960
		150875	Complete kit 480		480
		150877	Complete kit 96		96
		Apple stem pitting virus monoclonal - monoclonal	ASPV	151812	IgG
151815	IgG			0.1 ml	500
151819	IgG			0.025 ml	100
151822	Conjugate			0.2 ml	1000
151825	Conjugate			0.1 ml	500
151829	Conjugate			0.025 ml	100
151853	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
151862	Reagent set 960				960
151865	Reagent set 480				480
151872	Complete kit 960				960
151875	Complete kit 480				480
151877	Complete kit 96				96
Arabis mosaic virus polyclonal - polyclonal	ArMV			120312	IgG
		120315	IgG	0.1 ml	500
		120319	IgG	0.025 ml	100
		120322	Conjugate	0.2 ml	1000
		120325	Conjugate	0.1 ml	500
		120329	Conjugate	0.025 ml	100
		120353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		120362	Reagent set 960		960
		120365	Reagent set 480		480
		120372	Complete kit 960		960
		120373	Grapevine Complete kit 960		960
		120375	Complete kit 480		480
		120376	Grapevine Complete kit 480		480
		120377	Complete kit 96		96
		120378	Grapevine Complete kit 96		96
Arabis mosaic virus + Grapevine fanleaf virus polyclonal / monoclonal - monoclonal	ArMV + GFLV	120212	IgG	0.2 ml	1000
		120215	IgG	0.1 ml	500
		120219	IgG	0.025 ml	100
		120222	Conjugate	0.2 ml	1000
		120225	Conjugate	0.1 ml	500
		120229	Conjugate	0.025 ml	100
		120253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		120262	Reagent set 960		960
		120265	Reagent set 480		480
		120273	Complete kit 960		960
		120276	Complete kit 480		480
120278	Complete kit 96		96		

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Barley yellow dwarf virus - B (PAV-like) polyclonal - polyclonal	BYDV-B	140112	IgG	0.2 ml	1000
		140115	IgG	0.1 ml	500
		140119	IgG	0.025 ml	100
		140122	Conjugate	0.2 ml	1000
		140125	Conjugate	0.1 ml	500
		140129	Conjugate	0.025 ml	100
		140153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		140162	Reagent set 960		960
		140165	Reagent set 480		480
		140172	Complete kit 960		960
		140175	Complete kit 480		480
		140177	Complete kit 96		96
		Barley yellow dwarf virus - F (MAV-like) polyclonal - polyclonal	BYDV - F	140212	IgG
140215	IgG			0.1 ml	500
140219	IgG			0.025 ml	100
140222	Conjugate			0.2 ml	1000
140225	Conjugate			0.1 ml	500
140229	Conjugate			0.025 ml	100
140253	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
140262	Reagent set 960				960
140265	Reagent set 480				480
140272	Complete kit 960				960
140275	Complete kit 480				480
140277	Complete kit 96				96
Barley yellow dwarf virus - RPV polyclonal - polyclonal Synonym: Cereal yellow dwarf virus - RPV (CYDV - RPV)	BYDV - RPV			140612	IgG
		140615	IgG	0.1 ml	500
		140619	IgG	0.025 ml	100
		140622	Conjugate	0.2 ml	1000
		140625	Conjugate	0.1 ml	500
		140629	Conjugate	0.025 ml	100
		140653	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		140662	Reagent set 960		960
		140665	Reagent set 480		480
		140672	Complete kit 960		960
		140675	Complete kit 480		480
		140677	Complete kit 96		96
		Bean common mosaic virus polyclonal - polyclonal	BCMV	162012	IgG
162015	IgG			0.1 ml	500
162019	IgG			0.025 ml	100
162022	Conjugate			0.2 ml	1000
162025	Conjugate			0.1 ml	500
162029	Conjugate			0.025 ml	100
162053	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
162062	Reagent set 960				960
162065	Reagent set 480				480
162072	Complete kit 960				960
162075	Complete kit 480				480
162077	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Bean common mosaic necrosis virus polyclonal - polyclonal	BCMNV	162112	IgG	0.2 ml	1000
		162115	IgG	0.1 ml	500
		162119	IgG	0.025 ml	100
		162122	Conjugate	0.2 ml	1000
		162125	Conjugate	0.1 ml	500
		162129	Conjugate	0.025 ml	100
		162153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		162162	Reagent set 960		960
		162165	Reagent set 480		480
		162172	Complete kit 960		960
		162175	Complete kit 480		480
		162177	Complete kit 96		96
		Bean golden mosaic virus	BGMV	see Tomato yellow leaf curl virus	
Beet necrotic yellow vein virus (Rhizomania) polyclonal - polyclonal	BNYVV	160112	IgG	0.2 ml	1000
		160115	IgG	0.1 ml	500
		160119	IgG	0.025 ml	100
		160122	Conjugate	0.2 ml	1000
		160125	Conjugate	0.1 ml	500
		160129	Conjugate	0.025 ml	100
		160153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160162	Reagent set 960		960
		160165	Reagent set 480		480
		160172	Complete kit 960		960
		160175	Complete kit 480		480
		160177	Complete kit 96		96
		Calibrachoa mottle virus polyclonal - polyclonal Please note: available on request	CbMV	170712	IgG
170715	IgG			0.1 ml	500
170719	IgG			0.025 ml	100
170722	Conjugate			0.2 ml	1000
170725	Conjugate			0.1 ml	500
170729	Conjugate			0.025 ml	100
170753	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170762	Reagent set 960				960
170765	Reagent set 480				480
170772	Complete kit 960				960
170775	Complete kit 480				480
170777	Complete kit 96				96
Cauliflower mosaic virus polyclonal - polyclonal Needs special extraction buffer «Caulimo»	CaMV			162212	IgG
		162215	IgG	0.1 ml	500
		162219	IgG	0.025 ml	100
		162222	Conjugate	0.2 ml	1000
		162225	Conjugate	0.1 ml	500
		162229	Conjugate	0.025 ml	100
		162253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		162262	Reagent set 960		960
		162265	Reagent set 480		480
		162272	Complete kit 960		960
		162275	Complete kit 480		480
		162277	Complete kit 96		96
		Cereal yellow dwarf virus - RPV	CYDV - RPV	see Barley yellow dwarf virus - RPV	

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Cherry leaf roll virus polyclonal - polyclonal Suitable for cherry, peach, walnut and birch isolates.	CLRV	150912	IgG	0.2 ml	1000
		150915	IgG	0.1 ml	500
		150919	IgG	0.025 ml	100
		150922	Conjugate	0.2 ml	1000
		150925	Conjugate	0.1 ml	500
		150929	Conjugate	0.025 ml	100
		150953	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		150962	Reagent set 960		960
		150965	Reagent set 480		480
		150972	Complete kit 960		960
		150975	Complete kit 480		480
		150977	Complete kit 96		96
		Citrus tristeza virus polyclonal / monoclonal - polyclonal / monoclonal	CTV	151512	IgG
151515	IgG			0.1 ml	500
151519	IgG			0.025 ml	100
151522	Conjugate			0.2 ml	1000
151525	Conjugate			0.1 ml	500
151529	Conjugate			0.025 ml	100
151553	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
151562	Reagent set 960				960
151565	Reagent set 480				480
151572	Complete kit 960				960
151575	Complete kit 480				480
151577	Complete kit 96				96
Cucumber green mottle mosaic virus polyclonal - polyclonal Please note: available on request	CGMMV			160712	IgG
		160715	IgG	0.1 ml	500
		160719	IgG	0.025 ml	100
		160722	Conjugate	0.2 ml	1000
		160725	Conjugate	0.1 ml	500
		160729	Conjugate	0.025 ml	100
		160753	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160762	Reagent set 960		960
		160765	Reagent set 480		480
		160772	Complete kit 960		960
		160775	Complete kit 480		480
		160777	Complete kit 96		96
		Cucumber mosaic virus polyclonal - polyclonal / monoclonal	CMV	160612	IgG
160615	IgG			0.1 ml	500
160619	IgG			0.025 ml	100
160622	Conjugate			0.2 ml	1000
160625	Conjugate			0.1 ml	500
160629	Conjugate			0.025 ml	100
160653	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
160662	Reagent set 960				960
160665	Reagent set 480				480
160672	Complete kit 960				960
160675	Complete kit 480				480
160677	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Cucurbit yellow stunting disorder virus polyclonal - polyclonal	CYSDV	162312	IgG	0.2 ml	1000
		162315	IgG	0.1 ml	500
		162319	IgG	0.025 ml	100
		162322	Conjugate	0.2 ml	1000
		162325	Conjugate	0.1 ml	500
		162329	Conjugate	0.025 ml	100
		162353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		162362	Reagent set 960		960
		162365	Reagent set 480		480
		162372	Complete kit 960		960
		162375	Complete kit 480		480
		162377	Complete kit 96		96
		Erysimum latent virus polyclonal - polyclonal	ErLV	170912	IgG
170915	IgG			0.1 ml	500
170919	IgG			0.025 ml	100
170922	Conjugate			0.2 ml	1000
170925	Conjugate			0.1 ml	500
170929	Conjugate			0.025 ml	100
170953	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170962	Reagent set 960				960
170965	Reagent set 480				480
170972	Complete kit 960				960
170975	Complete kit 480				480
170977	Complete kit 96				96
Garlic common latent virus polyclonal - polyclonal	GCLV			160412	IgG
		160415	IgG	0.1 ml	500
		160419	IgG	0.025 ml	100
		160422	Conjugate	0.2 ml	1000
		160425	Conjugate	0.1 ml	500
		160429	Conjugate	0.025 ml	100
		160453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160462	Reagent set 960		960
		160465	Reagent set 480		480
		160472	Complete kit 960		960
		160475	Complete kit 480		480
		160477	Complete kit 96		96
		Grapevine fanleaf virus polyclonal - monoclonal	GFLV	120412	IgG
120415	IgG			0.1 ml	500
120419	IgG			0.025 ml	100
120422	Conjugate			0.2 ml	1000
120425	Conjugate			0.1 ml	500
120429	Conjugate			0.025 ml	100
120453	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
120462	Reagent set 960				960
120465	Reagent set 480				480
120473	Complete kit 960				960
120476	Complete kit 480				480
120478	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Grapevine fleck virus polyclonal - monoclonal	GFkV	120612	IgG	0.2 ml	1000
		120615	IgG	0.1 ml	500
		120619	IgG	0.025 ml	100
		120622	Conjugate	0.2 ml	1000
		120625	Conjugate	0.1 ml	500
		120629	Conjugate	0.025 ml	100
		120653	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		120662	Reagent set 960		960
		120665	Reagent set 480		480
		120673	Complete kit 960		960
		120676	Complete kit 480		480
		120678	Complete kit 96		96
		Grapevine leafroll - assoc. virus 1 monoclonal - monoclonal	GLRaV-1	120112	IgG
120115	IgG			0.1 ml	500
120119	IgG			0.025 ml	100
120122	Conjugate			0.2 ml	1000
120125	Conjugate			0.1 ml	500
120129	Conjugate			0.025 ml	100
120153	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
120162	Reagent set 960				960
120165	Reagent set 480				480
120173	Complete kit 960				960
120176	Complete kit 480				480
120178	Complete kit 96				96
Grapevine leafroll - assoc. virus 1+3 polyclonal / monoclonal - monoclonal	GLRaV-1+3			120912	IgG
		120915	IgG	0.1 ml	500
		120919	IgG	0.025 ml	100
		120922	Conjugate	0.2 ml	1000
		120925	Conjugate	0.1 ml	500
		120929	Conjugate	0.025 ml	100
		120953	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		120962	Reagent set 960		960
		120965	Reagent set 480		480
		120973	Complete kit 960		960
		120976	Complete kit 480		480
		120978	Complete kit 96		96
		Grapevine leafroll - assoc. virus 2 monoclonal - polyclonal	GLRaV-2	120812	IgG
120815	IgG			0.1 ml	500
120819	IgG			0.025 ml	100
120822	Conjugate			0.2 ml	1000
120825	Conjugate			0.1 ml	500
120829	Conjugate			0.025 ml	100
120853	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
120862	Reagent set 960				960
120865	Reagent set 480				480
120873	Complete kit 960				960
120876	Complete kit 480				480
120878	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Grapevine leafroll - assoc. virus 3 polyclonal / monoclonal - monoclonal	GLRaV-3	122312	IgG	0.2 ml	1000
		122315	IgG	0.1 ml	500
		122319	IgG	0.025 ml	100
		122322	Conjugate	0.2 ml	1000
		122325	Conjugate	0.1 ml	500
		122329	Conjugate	0.025 ml	100
		122353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		122362	Reagent set 960		960
		122365	Reagent set 480		480
		122373	Complete kit 960		960
		122376	Complete kit 480		480
		122378	Complete kit 96		96
		Grapevine leafroll - assoc. virus generic 4 strains (former GLRaV-4-9) monoclonal - monoclonal Recognizes GLRaV-4, GLRaV-4 strain 5, GLRaV-4 strain 6, GLRaV-4 strain 9 and GLRaV-4 strain Ob	GLRaV-4 strains	122512	IgG
122515	IgG			0.1 ml	500
122519	IgG			0.025 ml	100
122522	Conjugate			0.2 ml	1000
122525	Conjugate			0.1 ml	500
122529	Conjugate			0.025 ml	100
122553	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
122562	Reagent set 960				960
122565	Reagent set 480				480
122573	Complete kit 960				960
122576	Complete kit 480				480
122578	Complete kit 96				96
Grapevine leafroll - assoc. virus 4 strain 6 (former GLRaV-6) polyclonal - monoclonal	GLRaV-4 strain 6			122412	IgG
		122415	IgG	0.1 ml	500
		122419	IgG	0.025 ml	100
		122422	Conjugate	0.2 ml	1000
		122425	Conjugate	0.1 ml	500
		122429	Conjugate	0.025 ml	100
		122453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		122462	Reagent set 960		960
		122465	Reagent set 480		480
		122473	Complete kit 960		960
		122476	Complete kit 480		480
		122478	Complete kit 96		96
		Grapevine pinot gris virus polyclonal - polyclonal	GPGV	120712	IgG
120715	IgG			0.1 ml	500
120719	IgG			0.025 ml	100
120722	Conjugate			0.2 ml	1000
120725	Conjugate			0.1 ml	500
120729	Conjugate			0.025 ml	100
120753	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
120762	Reagent set 960				960
120765	Reagent set 480				480
120773	Complete kit 960				960
120776	Complete kit 480				480
120778	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Grapevine virus A polyclonal - polyclonal / monoclonal	GVA	122212	IgG	0.2 ml	1000
		122215	IgG	0.1 ml	500
		122219	IgG	0.025 ml	100
		122222	Conjugate	0.2 ml	1000
		122225	Conjugate	0.1 ml	500
		122229	Conjugate	0.025 ml	100
		122253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		122262	Reagent set 960		960
		122265	Reagent set 480		480
		122273	Complete kit 960		960
		122276	Complete kit 480		480
		122278	Complete kit 96		96
Groundnut ringspot virus	GRSV	see Tomato spotted wilt virus			
Impatiens necrotic spot virus polyclonal - monoclonal	INSV	190212	IgG	0.2 ml	1000
		190215	IgG	0.1 ml	500
		190219	IgG	0.025 ml	100
		190222	Conjugate	0.2 ml	1000
		190225	Conjugate	0.1 ml	500
		190229	Conjugate	0.025 ml	100
		190253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		190262	Reagent set 960		960
		190265	Reagent set 480		480
		190272	Complete kit 960		960
		190275	Complete kit 480		480
		190277	Complete kit 96		96
Leek yellow stripe virus polyclonal - polyclonal	LYSV	160312	IgG	0.2 ml	1000
		160315	IgG	0.1 ml	500
		160319	IgG	0.025 ml	100
		160322	Conjugate	0.2 ml	1000
		160325	Conjugate	0.1 ml	500
		160329	Conjugate	0.025 ml	100
		160353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160362	Reagent set 960		960
		160365	Reagent set 480		480
		160372	Complete kit 960		960
		160375	Complete kit 480		480
		160377	Complete kit 96		96
Lettuce mosaic virus polyclonal - monoclonal For leaves and seeds	LMV	160512	IgG	0.2 ml	1000
		160515	IgG	0.1 ml	500
		160519	IgG	0.025 ml	100
		160522	Conjugate	0.2 ml	1000
		160525	Conjugate	0.1 ml	500
		160529	Conjugate	0.025 ml	100
		160553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160562	Reagent set 960		960
		160565	Reagent set 480		480
		160572	Complete kit 960		960
		160575	Complete kit 480		480
		160577	Complete kit 96		96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Maize chlorotic mottle virus polyclonal - polyclonal	MCMV	140712	IgG	0.2 ml	1000
		140715	IgG	0.1 ml	500
		140719	IgG	0.025 ml	100
		140722	Conjugate	0.2 ml	1000
		140725	Conjugate	0.1 ml	500
		140729	Conjugate	0.025 ml	100
		140753	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		140762	Reagent set 960		960
		140765	Reagent set 480		480
		140772	Complete kit 960		960
		140775	Complete kit 480		480
		140777	Complete kit 96		96
		Maize dwarf mosaic virus polyclonal - polyclonal	MDMV	140312	IgG
140315	IgG			0.1 ml	500
140319	IgG			0.025 ml	100
140322	Conjugate			0.2 ml	1000
140325	Conjugate			0.1 ml	500
140329	Conjugate			0.025 ml	100
140353	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
140362	Reagent set 960				960
140365	Reagent set 480				480
140372	Complete kit 960				960
140375	Complete kit 480				480
140377	Complete kit 96				96
Mirafiori lettuce big-vein virus polyclonal - polyclonal Please note: available on request	MiLBVV			161612	IgG
		161615	IgG	0.1 ml	500
		161619	IgG	0.025 ml	100
		161622	Conjugate	0.2 ml	1000
		161625	Conjugate	0.1 ml	500
		161629	Conjugate	0.025 ml	100
		161653	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161662	Reagent set 960		960
		161665	Reagent set 480		480
		161672	Complete kit 960		960
		161675	Complete kit 480		480
		161677	Complete kit 96		96
		Onion yellow dwarf virus polyclonal - polyclonal	OYDV	160212	IgG
160215	IgG			0.1 ml	500
160219	IgG			0.025 ml	100
160222	Conjugate			0.2 ml	1000
160225	Conjugate			0.1 ml	500
160229	Conjugate			0.025 ml	100
160253	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
160262	Reagent set 960				960
160265	Reagent set 480				480
160272	Complete kit 960				960
160275	Complete kit 480				480
160277	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Papaya ringspot virus (= Watermelon mosaic virus 1) polyclonal - polyclonal	PRSV	151912	IgG	0.2 ml	1000
		151915	IgG	0.1 ml	500
		151919	IgG	0.025 ml	100
		151922	Conjugate	0.2 ml	1000
		151925	Conjugate	0.1 ml	500
		151929	Conjugate	0.025 ml	100
		151953	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		151962	Reagent set 960		960
		151965	Reagent set 480		480
		151972	Complete kit 960		960
		151975	Complete kit 480		480
		151977	Complete kit 96		96
		Pelargonium flower break virus polyclonal - polyclonal Please note: available on request	PFBV	170212	IgG
170215	IgG			0.1 ml	500
170219	IgG			0.025 ml	100
170222	Conjugate			0.2 ml	1000
170225	Conjugate			0.1 ml	500
170229	Conjugate			0.025 ml	100
170253	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170262	Reagent set 960				960
170265	Reagent set 480				480
170272	Complete kit 960				960
170275	Complete kit 480				480
170277	Complete kit 96				96
Pelargonium leaf curl virus polyclonal - polyclonal Please note: available on request	PLCV			170312	IgG
		170315	IgG	0.1 ml	500
		170319	IgG	0.025 ml	100
		170322	Conjugate	0.2 ml	1000
		170325	Conjugate	0.1 ml	500
		170329	Conjugate	0.025 ml	100
		170353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		170362	Reagent set 960		960
		170365	Reagent set 480		480
		170372	Complete kit 960		960
		170375	Complete kit 480		480
		170377	Complete kit 96		96
		Pelargonium line pattern virus polyclonal - polyclonal Please note: available on request	PLPV	170412	IgG
170415	IgG			0.1 ml	500
170419	IgG			0.025 ml	100
170422	Conjugate			0.2 ml	1000
170425	Conjugate			0.1 ml	500
170429	Conjugate			0.025 ml	100
170453	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170462	Reagent set 960				960
170465	Reagent set 480				480
170472	Complete kit 960				960
170475	Complete kit 480				480
170477	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Pepino mosaic virus polyclonal - polyclonal	PepMV	161512	IgG	0.2 ml	1000
		161515	IgG	0.1 ml	500
		161519	IgG	0.025 ml	100
		161522	Conjugate	0.2 ml	1000
		161525	Conjugate	0.1 ml	500
		161529	Conjugate	0.025 ml	100
		161553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161562	Reagent set 960		960
		161565	Reagent set 480		480
		161572	Complete kit 960		960
		161575	Complete kit 480		480
		161577	Complete kit 96		96
		Pepper mild mottle virus polyclonal - polyclonal	PMMoV	161712	IgG
161715	IgG			0.1 ml	500
161719	IgG			0.025 ml	100
161722	Conjugate			0.2 ml	1000
161725	Conjugate			0.1 ml	500
161729	Conjugate			0.025 ml	100
161753	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
161762	Reagent set 960				960
161765	Reagent set 480				480
161772	Complete kit 960				960
161775	Complete kit 480				480
161777	Complete kit 96				96
Petunia asteroid mosaic virus polyclonal - polyclonal	PeAMV			170812	IgG
		170815	IgG	0.1 ml	500
		170819	IgG	0.025 ml	100
		170822	Conjugate	0.2 ml	1000
		170825	Conjugate	0.1 ml	500
		170829	Conjugate	0.025 ml	100
		170853	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		170862	Reagent set 960		960
		170865	Reagent set 480		480
		170872	Complete kit 960		960
		170875	Complete kit 480		480
		170877	Complete kit 96		96
		Plum pox virus (Sharka) polyclonal - polyclonal / monoclonal	PPV	150512	IgG
150515	IgG			0.1 ml	500
150519	IgG			0.025 ml	100
150522	Conjugate			0.2 ml	1000
150525	Conjugate			0.1 ml	500
150529	Conjugate			0.025 ml	100
150553	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
150562	Reagent set 960				960
150565	Reagent set 480				480
150572	Complete kit 960				960
150575	Complete kit 480				480
150577	Complete kit 96				96

Please note: available on request

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Potato leafroll virus polyclonal - monoclonal	PLRV	110611	IgG	1.0 ml	5000
		110612	IgG	0.2 ml	1000
		110615	IgG	0.1 ml	500
		110619	IgG	0.025 ml	100
		110621	Conjugate	1.0 ml	5000
		110622	Conjugate	0.2 ml	1000
		110625	Conjugate	0.1 ml	500
		110629	Conjugate	0.025 ml	100
		110653	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		110670	Reagent set 5000		5000
		110662	Reagent set 960		960
		110665	Reagent set 480		480
		110680	Complete kit 5000		5000
		110672	Complete kit 960		960
		110675	Complete kit 480		480
		110677	Complete kit 96		96
		Potato mop-top virus monoclonal - monoclonal	PMTV	113012	IgG
113015	IgG			0.1 ml	500
113019	IgG			0.025 ml	100
113022	Conjugate			0.2 ml	1000
113025	Conjugate			0.1 ml	500
113029	Conjugate			0.025 ml	100
113053	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
113062	Reagent set 960				960
113065	Reagent set 480				480
113072	Complete kit 960				960
113075	Complete kit 480				480
113077	Complete kit 96				96
Potato virus A monoclonal - monoclonal	PVA	112111	IgG	1.0 ml	5000
		112112	IgG	0.2 ml	1000
		112115	IgG	0.1 ml	500
		112119	IgG	0.025 ml	100
		112121	Conjugate	1.0 ml	5000
		112122	Conjugate	0.2 ml	1000
		112125	Conjugate	0.1 ml	500
		112129	Conjugate	0.025 ml	100
		112153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		112170	Reagent set 5000		5000
		112162	Reagent set 960		960
		112165	Reagent set 480		480
		112180	Complete kit 5000		5000
		112172	Complete kit 960		960
		112175	Complete kit 480		480
		112177	Complete kit 96		96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Potato virus M polyclonal - polyclonal	PVM	110211	IgG	1.0 ml	5000
		110212	IgG	0.2 ml	1000
		110215	IgG	0.1 ml	500
		110219	IgG	0.025 ml	100
		110221	Conjugate	1.0 ml	5000
		110222	Conjugate	0.2 ml	1000
		110225	Conjugate	0.1 ml	500
		110229	Conjugate	0.025 ml	100
		110253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		110270	Reagent set 5000		5000
		110262	Reagent set 960		960
		110265	Reagent set 480		480
		110280	Complete kit 5000		5000
		110272	Complete kit 960		960
		110275	Complete kit 480		480
		110277	Complete kit 96		96
		Potato virus S polyclonal - polyclonal	PVS	110311	IgG
110312	IgG			0.2 ml	1000
110315	IgG			0.1 ml	500
110319	IgG			0.025 ml	100
110321	Conjugate			1.0 ml	5000
110322	Conjugate			0.2 ml	1000
110325	Conjugate			0.1 ml	500
110329	Conjugate			0.025 ml	100
110353	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
110370	Reagent set 5000				5000
110362	Reagent set 960				960
110365	Reagent set 480				480
110380	Complete kit 5000				5000
110372	Complete kit 960				960
110375	Complete kit 480				480
110377	Complete kit 96				96
Potato virus V (PVY^c 624) polyclonal - polyclonal	PVV			110812	IgG
		110815	IgG	0.1 ml	500
		110819	IgG	0.025 ml	100
		110822	Conjugate	0.2 ml	1000
		110825	Conjugate	0.1 ml	500
		110829	Conjugate	0.025 ml	100
		110853	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		110862	Reagent set 960		960
		110865	Reagent set 480		480
		110872	Complete kit 960		960
		110875	Complete kit 480		480
110877	Complete kit 96		96		

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Potato virus X polyclonal - polyclonal	PVX	110411	IgG	1.0 ml	5000
		110412	IgG	0.2 ml	1000
		110415	IgG	0.1 ml	500
		110419	IgG	0.025 ml	100
		110421	Conjugate	1.0 ml	5000
		110422	Conjugate	0.2 ml	1000
		110425	Conjugate	0.1 ml	500
		110429	Conjugate	0.025 ml	100
		110453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		110470	Reagent set 5000		5000
		110462	Reagent set 960		960
		110465	Reagent set 480		480
		110480	Complete kit 5000		5000
		110472	Complete kit 960		960
		110475	Complete kit 480		480
		110477	Complete kit 96		96
		Potato virus Y - monoclonal* monoclonal - monoclonal Traditional broad-spectrum reagent. Recognizes all isolates except atypical isolates of PVY ^{bb} subgroup *For routine diagnostics, we recommend our full-spectrum PVY reagent ("PVY - monoclonal cocktail").	PVY-mono	112511	IgG
112512	IgG			0.2 ml	1000
112515	IgG			0.1 ml	500
112519	IgG			0.025 ml	100
112521	Conjugate			1.0 ml	5000
112522	Conjugate			0.2 ml	1000
112525	Conjugate			0.1 ml	500
112529	Conjugate			0.025 ml	100
112553	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
112570	Reagent set 5000				5000
112562	Reagent set 960				960
112565	Reagent set 480				480
112580	Complete kit 5000				5000
112572	Complete kit 960				960
112575	Complete kit 480				480
112577	Complete kit 96				96
Potato virus Y - monoclonal cocktail monoclonal - monoclonal Full-spectrum reagent. Contains complementary monoclonal antibodies to isolates from different strain groups. Recognizes all known isolates of all groups/subgroups of PVY."	PVY- monocock			112911	IgG
		112912	IgG	0.2 ml	1000
		112915	IgG	0.1 ml	500
		112919	IgG	0.025 ml	100
		112921	Conjugate	1.0 ml	5000
		112922	Conjugate	0.2 ml	1000
		112925	Conjugate	0.1 ml	500
		112929	Conjugate	0.025 ml	100
		112553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		112970	Reagent Set 5000		5000
		112962	Reagent set 960		960
		112965	Reagent set 480		480
		112980	Complete kit 5000		5000
		112972	Complete kit 960		960
		112975	Complete kit 480		480
		112977	Complete kit 96		96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays		
Potato virus Y - necrotic* monoclonal - monoclonal Specifically recognizes PVY necrotic strain group, except Wilga type (PVY ^N W). DAS-ELISA format is essential to ensure described N-specificity. *For routine diagnostics, we recommend our full-spectrum PVY reagent ("PVY - monoclonal cocktail").	PVY ^N	112712	IgG	0.2 ml	1000		
		112715	IgG	0.1 ml	500		
		112719	IgG	0.025 ml	100		
		112722	Conjugate	0.2 ml	1000		
		112725	Conjugate	0.1 ml	500		
		112729	Conjugate	0.025 ml	100		
		112553	Positive control	2.5 ml	12		
		180043	Negative control	2.5 ml	12		
		112762	Reagent set 960		960		
		112765	Reagent set 480		480		
		112772	Complete kit 960		960		
		112775	Complete kit 480		480		
		112777	Complete kit 96		96		
		Potato virus Y - polyclonal* polyclonal - polyclonal Serological spectrum is similar to PVY monoclonal cocktail. *For routine diagnostics, we recommend our full-spectrum PVY reagent ("PVY - monoclonal cocktail").	PVY - poly	110511	IgG	1.0 ml	5000
				110512	IgG	0.2 ml	1000
110515	IgG			0.1 ml	500		
110519	IgG			0.025 ml	100		
110521	Conjugate			1.0 ml	5000		
110522	Conjugate			0.2 ml	1000		
110525	Conjugate			0.1 ml	500		
110529	Conjugate			0.025 ml	100		
112553	Positive control			2.5 ml	12		
180043	Negative control			2.5 ml	12		
110570	Reagent set 5000				5000		
110562	Reagent set 960				960		
110565	Reagent set 480				480		
110580	Complete kit 5000				5000		
110572	Complete kit 960				960		
110575	Complete kit 480		480				
110577	Complete kit 96		96				
Poty group test monoclonal Only available as kit in indirect PTA-ELISA format.	Poty group	163072	Complete kit 960		960		
		163075	Complete kit 480		480		
Prune dwarf virus polyclonal - polyclonal	PDV	150612	IgG	0.2 ml	1000		
		150615	IgG	0.1 ml	500		
		150619	IgG	0.025 ml	100		
		150622	Conjugate	0.2 ml	1000		
		150625	Conjugate	0.1 ml	500		
		150629	Conjugate	0.025 ml	100		
		150653	Positive control	2.5 ml	12		
		180043	Negative control	2.5 ml	12		
		150662	Reagent set 960		960		
		150665	Reagent set 480		480		
		150672	Complete kit 960		960		
		150675	Complete kit 480		480		
150677	Complete kit 96		96				

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Prunus necrotic ringspot virus polyclonal - polyclonal	PNRSV	151112	IgG	0.2 ml	1000
		151115	IgG	0.1 ml	500
		151119	IgG	0.025 ml	100
		151122	Conjugate	0.2 ml	1000
		151125	Conjugate	0.1 ml	500
		151129	Conjugate	0.025 ml	100
		151153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		151162	Reagent set 960		960
		151165	Reagent set 480		480
		151172	Complete kit 960		960
		151175	Complete kit 480		480
		151177	Complete kit 96		96
		Radish mosaic virus polyclonal - polyclonal Please note: available on request	RaMV	162412	IgG
162415	IgG			0.1 ml	500
162419	IgG			0.025 ml	100
162422	Conjugate			0.2 ml	1000
162425	Conjugate			0.1 ml	500
162429	Conjugate			0.025 ml	100
162453	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
162462	Reagent set 960				960
162465	Reagent set 480				480
162472	Complete kit 960				960
162475	Complete kit 480				480
162477	Complete kit 96				96
Ralstonia solanacearum polyclonal - polyclonal	Rs			170512	IgG
		170515	IgG	0.1 ml	500
		170519	IgG	0.025 ml	100
		170522	Conjugate	0.2 ml	1000
		170525	Conjugate	0.1 ml	500
		170529	Conjugate	0.025 ml	100
		170553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		170562	Reagent set 960		960
		170565	Reagent set 480		480
		170572	Complete kit 960		960
		170575	Complete kit 480		480
		170577	Complete kit 96		96
		Raspberry bushy dwarf virus polyclonal - monoclonal	RBDV	152012	IgG
152015	IgG			0.1 ml	500
152019	IgG			0.025 ml	100
152022	Conjugate			0.2 ml	1000
152025	Conjugate			0.1 ml	500
152029	Conjugate			0.025 ml	100
152053	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
152062	Reagent set 960				960
152065	Reagent set 480				480
152072	Complete kit 960				960
152075	Complete kit 480				480
152077	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Raspberry ringspot virus - ch polyclonal - polyclonal Broad-spectrum reagent recognizes most common isolates from fruit tree, small fruit and grapevine.	RpRSV - ch	150412	IgG	0.2 ml	1000
		150415	IgG	0.1 ml	500
		150419	IgG	0.025 ml	100
		150422	Conjugate	0.2 ml	1000
		150425	Conjugate	0.1 ml	500
		150429	Conjugate	0.025 ml	100
		150453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		150462	Reagent set 960		960
		150465	Reagent set 480		480
		150472	Complete kit 960		960
		150473	Grapevine Complete kit 960		960
		150475	Complete kit 480		480
		150476	Grapevine Complete kit 480		480
		150477	Complete kit 96		96
			150478	Grapevine Complete kit 96	
Raspberry ringspot virus - g* polyclonal - polyclonal Specific reagent for certain grapevine isolates from limited geographic locations. Complementary to RpRSV-ch. *For routine diagnostics, we recommend our RpRSV-ch reagents.	RpRSV - g	120512	IgG	0.2 ml	1000
		120515	IgG	0.1 ml	500
		120519	IgG	0.025 ml	100
		120522	Conjugate	0.2 ml	1000
		120525	Conjugate	0.1 ml	500
		120529	Conjugate	0.025 ml	100
		120553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		120562	Reagent set 960		960
		120565	Reagent set 480		480
		120573	Complete kit 960		960
		120576	Complete kit 480		480
			120578	Complete kit 96	
Soybean mosaic virus polyclonal - polyclonal Please note: available on request	SMV	162612	IgG	0.2 ml	1000
		162615	IgG	0.1 ml	500
		162619	IgG	0.025 ml	100
		162622	Conjugate	0.2 ml	1000
		162625	Conjugate	0.1 ml	500
		162629	Conjugate	0.025 ml	100
		162653	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		162662	Reagent set 960		960
		162665	Reagent set 480		480
		162672	Complete kit 960		960
		162675	Complete kit 480		480
			162677	Complete kit 96	
Spongospora subterranea f. sp. subterranea (Powdery scab of potato) monoclonal - monoclonal Suitable for detecting resting spores on potato tubers and in soil.	Sss	111112	IgG	0.2 ml	1000
		111115	IgG	0.1 ml	500
		111119	IgG	0.025 ml	100
		111122	Conjugate	0.2 ml	1000
		111125	Conjugate	0.1 ml	500
		111129	Conjugate	0.025 ml	100
		111153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		111162	Reagent set 960		960
		111165	Reagent set 480		480
		111172	Complete kit 960		960
		111175	Complete kit 480		480
			111177	Complete kit 96	
Squash leaf curl virus	SLCV	see Tomato yellow leaf curl virus			

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Squash mosaic virus polyclonal - polyclonal	SqMV	160912	IgG	0.2 ml	1000
		160915	IgG	0.1 ml	500
		160919	IgG	0.025 ml	100
		160922	Conjugate	0.2 ml	1000
		160925	Conjugate	0.1 ml	500
		160929	Conjugate	0.025 ml	100
		160953	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		160962	Reagent set 960		960
		160965	Reagent set 480		480
		160972	Complete kit 960		960
		160975	Complete kit 480		480
		160977	Complete kit 96		96
		Strawberry latent ringspot virus polyclonal - polyclonal	SLRSV	151612	IgG
151615	IgG			0.1 ml	500
151619	IgG			0.025 ml	100
151622	Conjugate			0.2 ml	1000
151625	Conjugate			0.1 ml	500
151629	Conjugate			0.025 ml	100
151653	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
151662	Reagent set 960				960
151665	Reagent set 480				480
151672	Complete kit 960				960
151673	Grapevine Complete kit 960				960
151675	Complete kit 480				480
151676	Grapevine Complete kit 480				480
151677	Complete kit 96				96
Strawberry mild yellow edge potexvirus polyclonal - monoclonal	SMYEPV			152112	IgG
		152115	IgG	0.1 ml	500
		152119	IgG	0.025 ml	100
		152122	Conjugate	0.2 ml	1000
		152125	Conjugate	0.1 ml	500
		152129	Conjugate	0.025 ml	100
		152153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		152162	Reagent set 960		960
		152165	Reagent set 480		480
		152172	Complete kit 960		960
		152175	Complete kit 480		480
		152177	Complete kit 96		96
		Sugarcane mosaic virus polyclonal - polyclonal	SCMV	140412	IgG
140415	IgG			0.1 ml	500
140419	IgG			0.025 ml	100
140422	Conjugate			0.2 ml	1000
140425	Conjugate			0.1 ml	500
140429	Conjugate			0.025 ml	100
140453	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
140462	Reagent set 960				960
140465	Reagent set 480				480
140472	Complete kit 960				960
140475	Complete kit 480				480
140477	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Tobacco mosaic virus polyclonal - polyclonal	TMV	190412	IgG	0.2 ml	1000
		190415	IgG	0.1 ml	500
		190419	IgG	0.025 ml	100
		190422	Conjugate	0.2 ml	1000
		190425	Conjugate	0.1 ml	500
		190429	Conjugate	0.025 ml	100
		190453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		190462	Reagent set 960		960
		190465	Reagent set 480		480
		190472	Complete kit 960		960
		190475	Complete kit 480		480
		190477	Complete kit 96		96
		Tobacco necrosis virus polyclonal - polyclonal Please note: available on request	TNV	152412	IgG
152415	IgG			0.1 ml	500
152419	IgG			0.025 ml	100
152422	Conjugate			0.2 ml	1000
152425	Conjugate			0.1 ml	500
152429	Conjugate			0.025 ml	100
152453	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
152462	Reagent set 960				960
152465	Reagent set 480				480
152472	Complete kit 960				960
152475	Complete kit 480				480
152477	Complete kit 96				96
Tobacco ringspot virus polyclonal - polyclonal	TRSV			152212	IgG
		152215	IgG	0.1 ml	500
		152219	IgG	0.025 ml	100
		152222	Conjugate	0.2 ml	1000
		152225	Conjugate	0.1 ml	500
		152229	Conjugate	0.025 ml	100
		152253	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		152262	Reagent set 960		960
		152265	Reagent set 480		480
		152272	Complete kit 960		960
		152273	Grapevine Complete kit 960		960
		152275	Complete kit 480		480
		152276	Grapevine Complete kit 480		480
		152277	Complete kit 96		96
		152278	Grapevine Complete kit 96		96
Tobacco streak virus polyclonal - monoclonal	TSV	152312	IgG	0.2 ml	1000
		152315	IgG	0.1 ml	500
		152319	IgG	0.025 ml	100
		152322	Conjugate	0.2 ml	1000
		152325	Conjugate	0.1 ml	500
		152329	Conjugate	0.025 ml	100
		152353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		152362	Reagent set 960		960
		152365	Reagent set 480		480
		152372	Complete kit 960		960
		152375	Complete kit 480		480
		152377	Complete kit 96		96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Tomato apex necrosis virus Synonym: Tomato marchitez virus (ToMarV) polyclonal - polyclonal	ToANV	152712	IgG	0.2 ml	1000
		152715	IgG	0.1 ml	500
		152719	IgG	0.025 ml	100
		152722	Conjugate	0.2 ml	1000
		152725	Conjugate	0.1 ml	500
		152729	Conjugate	0.025 ml	100
		152753	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		152762	Reagent set 960		960
		152765	Reagent set 480		480
		152772	Complete kit 960		960
		152775	Complete kit 480		480
		152777	Complete kit 96		96
		Tomato black ring virus polyclonal / monoclonal - polyclonal / monoclonal Broad-spectrum reagent	TBRV	113112	IgG
113115	IgG			0.1 ml	500
113119	IgG			0.025 ml	100
113122	Conjugate			0.2 ml	1000
113125	Conjugate			0.1 ml	500
113129	Conjugate			0.025 ml	100
113153	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
113162	Reagent set 960				960
113165	Reagent set 480				480
113172	Complete kit 960				960
113173	Grapevine Complete kit 960				960
113175	Complete kit 480				480
113176	Grapevine Complete kit 480				480
113177	Complete kit 96				96
ToBRFV Tomato brown rugose fruit virus polyclonal - polyclonal Because of ToBRFV quarantine recommendations, the positive controls are made from extracts of <i>Tobamovirus</i> -infected plants (ToBRFV-negative).	ToBRFV			152812	IgG
		152815	IgG	0.1 ml	500
		152819	IgG	0.025 ml	100
		152822	Conjugate	0.2 ml	1000
		152825	Conjugate	0.1 ml	500
		152829	Conjugate	0.025 ml	100
		152853	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		152862	Reagent set 960		960
		152865	Reagent set 480		480
		152872	Complete kit 960		960
		152875	Complete kit 480		480
		152877	Complete kit 96		96
Tomato bushy stunt virus polyclonal - polyclonal Please note: available on request	TBSV	161812	IgG	0.2 ml	1000
		161815	IgG	0.1 ml	500
		161819	IgG	0.025 ml	100
		161822	Conjugate	0.2 ml	1000
		161825	Conjugate	0.1 ml	500
		161829	Conjugate	0.025 ml	100
		161853	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161862	Reagent set 960		960
		161865	Reagent set 480		480
		161872	Complete kit 960		960
161875	Complete kit 480		480		
161877	Complete kit 96		96		

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays	
Tomato chlorotic spot virus	TCSV	see Tomato spotted wilt virus				
Tomato mosaic virus	ToMV	152612	IgG	0.2 ml	1000	
polyclonal - polyclonal		152615	IgG	0.1 ml	500	
		152619	IgG	0.025 ml	100	
		152622	Conjugate	0.2 ml	1000	
		152625	Conjugate	0.1 ml	500	
		152629	Conjugate	0.025 ml	100	
		152653	Positive control	2.5 ml	12	
		180043	Negative control	2.5 ml	12	
		152662	Reagent set 960		960	
		152665	Reagent set 480		480	
		152672	Complete kit 960		960	
		152675	Complete kit 480		480	
		152677	Complete kit 96		96	
	Tomato mottle virus	ToMoV	see Tomato yellow leaf curl virus			
	Tomato ringspot virus	ToRSV	151412	IgG	0.2 ml	1000
polyclonal - polyclonal		151415	IgG	0.1 ml	500	
		151419	IgG	0.025 ml	100	
		151422	Conjugate	0.2 ml	1000	
		151425	Conjugate	0.1 ml	500	
		151429	Conjugate	0.025 ml	100	
		151453	Positive control	2.5 ml	12	
		180043	Negative control	2.5 ml	12	
		151462	Reagent set 960		960	
		151465	Reagent set 480		480	
		151472	Complete kit 960		960	
		151473	Grapevine Complete kit 960		960	
		151475	Complete kit 480		480	
		151476	Grapevine Complete kit 480		480	
		151477	Complete kit 96		96	
	151478	Grapevine Complete kit 96		96		
Tomato ringspot virus-Ch*	ToRSV-Ch	151312	IgG	0.2 ml	1000	
polyclonal - polyclonal		151315	IgG	0.1 ml	500	
		151319	IgG	0.025 ml	100	
		151322	Conjugate	0.2 ml	1000	
		151325	Conjugate	0.1 ml	500	
		151329	Conjugate	0.025 ml	100	
		151353	Positive control	2.5 ml	12	
		180043	Negative control	2.5 ml	12	
		151362	Reagent set 960		960	
		151365	Reagent set 480		480	
		151372	Complete kit 960		960	
		151373	Grapevine Complete kit 960		960	
		151375	Complete kit 480		480	
		151376	Grapevine Complete kit 480		480	
		151377	Complete kit 96		96	
	151378	Grapevine Complete kit 96		96		

Broad-spectrum reagent, recognizes most common isolates from many host plants (fruit trees, small fruits, grapevines, ornamentals).

Specific for Chickadee (originally from apple) and grape yellow vein (GYV) strains. Complementary to ToRSV-PYBM.

*For routine diagnostics, we recommend our ToRSV reagents.

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Tomato spotted wilt virus	TSWV	190112	IgG	0.2 ml	1000
		190115	IgG	0.1 ml	500
		190119	IgG	0.025 ml	100
		190122	Conjugate	0.2 ml	1000
		190125	Conjugate	0.1 ml	500
		190129	Conjugate	0.025 ml	100
		190153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		190162	Reagent set 960		960
		190165	Reagent set 480		480
		190172	Complete kit 960		960
		190175	Complete kit 480		480
		190177	Complete kit 96		96
polyclonal - monoclonal					
Recognizes serogroups I, II of Tospoviruses (TSWV, TCSV, GRSV).					
Tomato yellow leaf curl virus	TYLCV	161412	IgG	0.2 ml	1000
		161415	IgG	0.1 ml	500
		161419	IgG	0.025 ml	100
		161422	Conjugate	0.2 ml	1000
		161425	Conjugate	0.1 ml	500
		161429	Conjugate	0.025 ml	100
		161453	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161462	Reagent set 960		960
		161465	Reagent set 480		480
		161472	Complete kit 960		960
		161475	Complete kit 480		480
		161477	Complete kit 96		96
polyclonal - monoclonal					
Broad-spectrum, recognizes TYLCV (different isolates), African cassava mosaic virus (ACMV), Bean golden mosaic virus (BGMV), East African cassava mosaic virus EACMV), Indian cassava mosaic virus (ICMV), Squash leaf curl virus (SLCV) and Tomato mottle virus (ToMoV).					
Tospovirus broad-spectrum (serogroups I,II,III)	Tospo (I, II, III)	190312	IgG	0.2 ml	1000
		190315	IgG	0.1 ml	500
		190319	IgG	0.025 ml	100
		190322	Conjugate	0.2 ml	1000
		190325	Conjugate	0.1 ml	500
		190329	Conjugate	0.025 ml	100
		190353	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		190362	Reagent set 960		960
		190365	Reagent set 480		480
		190372	Complete kit 960		960
		190375	Complete kit 480		480
		190377	Complete kit 96		96
polyclonal - monoclonal					
Recognizes serogroups I, II,III of Tospoviruses (TSWV, TCSV, GRSV and INSV).					
Turnip mosaic virus	TuMV	161012	IgG	0.2 ml	1000
		161015	IgG	0.1 ml	500
		161019	IgG	0.025 ml	100
		161022	Conjugate	0.2 ml	1000
		161025	Conjugate	0.1 ml	500
		161029	Conjugate	0.025 ml	100
		161053	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161062	Reagent set 960		960
		161065	Reagent set 480		480
		161072	Complete kit 960		960
		161075	Complete kit 480		480
		161077	Complete kit 96		96
polyclonal - monoclonal					
Please note: available on request					

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays
Turnip yellow mosaic virus polyclonal - polyclonal	TYMV	162512	IgG	0.2 ml	1000
		162515	IgG	0.1 ml	500
		162519	IgG	0.025 ml	100
		162522	Conjugate	0.2 ml	1000
		162525	Conjugate	0.1 ml	500
		162529	Conjugate	0.025 ml	100
		162553	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		162562	Reagent set 960		960
		162565	Reagent set 480		480
		162572	Complete kit 960		960
		162575	Complete kit 480		480
		162577	Complete kit 96		96
		Verticillium spp. monoclonal - monoclonal	Verticillium	161312	IgG
161315	IgG			0.1 ml	500
161319	IgG			0.025 ml	100
161322	Conjugate			0.2 ml	1000
161325	Conjugate			0.1 ml	500
161329	Conjugate			0.025 ml	100
161353	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
161362	Reagent set 960				960
161365	Reagent set 480				480
161372	Complete kit 960				960
161375	Complete kit 480				480
161377	Complete kit 96				96
Watermelon mosaic virus 1	WMV - 1			see Papaya ringspot virus	
Watermelon mosaic virus 2 polyclonal - polyclonal	WMV - 2	161112	IgG	0.2 ml	1000
		161115	IgG	0.1 ml	500
		161119	IgG	0.025 ml	100
		161122	Conjugate	0.2 ml	1000
		161125	Conjugate	0.1 ml	500
		161129	Conjugate	0.025 ml	100
		161153	Positive control	2.5 ml	12
		180043	Negative control	2.5 ml	12
		161162	Reagent set 960		960
		161165	Reagent set 480		480
		161172	Complete kit 960		960
		161175	Complete kit 480		480
		161177	Complete kit 96		96
		Xanthomonas campestris pv. pelargonii Synonym: Xanthomonas hortorum (Xh) polyclonal - polyclonal	Xcp	170112	IgG
170115	IgG			0.1 ml	500
170119	IgG			0.025 ml	100
170122	Conjugate			0.2 ml	1000
170125	Conjugate			0.1 ml	500
170129	Conjugate			0.025 ml	100
170153	Positive control			2.5 ml	12
180043	Negative control			2.5 ml	12
170162	Reagent set 960				960
170165	Reagent set 480				480
170172	Complete kit 960				960
170175	Complete kit 480				480
170177	Complete kit 96				96

Pathogen	Abbr.	Art. No.	Reagent form	Size	Assays	
<i>Xanthomonas fragariae</i> polyclonal - polyclonal	Xfr	150312	IgG	0.2 ml	1000	
		150315	IgG	0.1 ml	500	
		150319	IgG	0.025 ml	100	
		150322	Conjugate	0.2 ml	1000	
		150325	Conjugate	0.1 ml	500	
		150329	Conjugate	0.025 ml	100	
		150353	Positive control	2.5 ml	12	
		180043	Negative control	2.5 ml	12	
		150362	Reagent set 960		960	
		150365	Reagent set 480		480	
		150372	Complete kit 960		960	
		150375	Complete kit 480		480	
		150377	Complete kit 96		96	
		<i>Xanthomonas hortorum</i>	Xh	see <i>Xanthomonas campestris</i> pv. <i>pelargonii</i>		
		Zucchini yellow mosaic virus polyclonal - polyclonal	ZYMV	161212	IgG	0.2 ml
161215	IgG			0.1 ml	500	
161219	IgG			0.025 ml	100	
161222	Conjugate			0.2 ml	1000	
161225	Conjugate			0.1 ml	500	
161229	Conjugate			0.025 ml	100	
161253	Positive control			2.5 ml	12	
180043	Negative control			2.5 ml	12	
161262	Reagent set 960				960	
161265	Reagent set 480				480	
161272	Complete kit 960				960	
161275	Complete kit 480				480	
161277	Complete kit 96				96	

qPCR Sets and Kits

By using BIOREBA's ready-to-use qPCR tests, customers **save** days of work compared to doing it themselves in the laboratory. Our qPCR products are supplied with lot specific data sheets. We are committed to **lot-to-lot consistency** for increased test reliability. Our qPCR sets include: Taq Master Mix, RT Master Mix, nuclease-free water, primers, probes and internal control mix.

qPCR PLRV/PVY and qPCR PVA/PVM

qPCR set and kit for the detection of Potato Viruses



The qPCR sets enable the simultaneous detection (multiplex) of potato virus Y (**PVY**; all strains) and potato leafroll virus (**PLRV**) or **PVA** (potato virus A) and **PVM** (potato virus M) respectively with real time PCR.

Suitable tissues for the source of viral RNA are potato leaves and tubers (peel, stolon, sprouts). It is highly recommended to use the BIOREBA "Potato DNA/RNA rapid extraction" (included in the qPCR kits) for the extraction of RNA. The qPCR for the detection of PLRV/PVY/IC and PVA/PVM/IC was developed and validated in collaboration with Agroscope, the Swiss centre of excellence for research in the agriculture and food sector.

There are two set sizes available for 96 and for 192 samples and there are four kit sizes available for pools of 10 or up to 25 potato samples and for 96 or 192 extractions respectively.

qPCR GRBV



qPCR set for the detection of Grapevine Red Blotch Virus (GRBV)

The qPCR set enables the simultaneous detection (multiplex) of Grapevine red blotch virus (**GRBV**) and a *Vitis vinifera* host gene control (IC) with real-time PCR.

Suitable tissue sources of viral DNA are grapevine leaves or wood scrapings (phloem). BIOREBA recommends to homogenize samples using the BIOREBA homogenizer system HOMEX and BIOREBA extraction bags Universal, and to follow a CTAB buffer-based method for DNA extraction.

The qPCR for the detection of GRBV was developed in collaboration with Agroscope, the Swiss center of excellence for research in the agriculture and food sectors.

There are two set sizes available for 96 and 192 samples.

qPCR ToBRFV



qPCR set for the detection of Tomato Brown Rugose Fruit Virus (ToBRFV)

The qPCR set enables the detection of Tomato brown rugose fruit virus (**ToBRFV**) with a triplex real-time RT-PCR.

ToBRFV with its resistance-breaking ability is responsible for enormous economic losses. The EPPO (European and Mediterranean Plant Protection Organization) recommends ToBRFV to be regulated as a quarantine pathogen. ToBRFV causes different symptoms in different varieties. Observed symptoms include yellow or brown, wrinkled spots on fruits, narrowing and mosaic discolorations of leaves as well as deformed or irregularly ripened fruits. This test is suitable for detecting ToBRFV in plant leaves, seeds and fruits.

There are two set sizes available for 96 and 192 samples.

Art. No.	Product	No. of reactions	Amount of samples (potatoes)
839600	qPCR PLRV/PVY set 96 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes/IC Mix_PLRV/PVY • Nuclease free water • Positive and negative controls 	96	
839200	qPCR PLRV/PVY set 192 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes/IC Mix_PLRV/PVY • Nuclease free water • Positive and negative controls 	192	
839610	qPCR PLRV/PVY kit 96/10 includes: <ul style="list-style-type: none"> • qPCR PLRV/PVY set 96 • Potato DNA/RNA rapid extraction set 10x for 96 extractions and a pool of 10 potatoes 	96	10
839625	qPCR PLRV/PVY kit 96/25 includes: <ul style="list-style-type: none"> • qPCR PLRV/PVY set 96 • Potato DNA/RNA rapid extraction set 25x for 96 extractions and a pool of 25 potatoes 	96	25
839210	qPCR PLRV/PVY kit 192/10 includes: <ul style="list-style-type: none"> • qPCR PLRV/PVY set 192 • Potato DNA/RNA rapid extraction set 10x for 192 extractions and a pool of 10 potatoes 	192	10
839225	qPCR PLRV/PVY kit 192/25 includes: <ul style="list-style-type: none"> • qPCR PLRV/PVY set 192 • Potato DNA/RNA rapid extraction set 25x for 192 extractions and a pool of 25 potatoes 	192	25
849600	qPCR PVA/PVM set 96 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes/IC Mix_PVA/PVM/IC • Nuclease free water • Positive and negative controls 	96	
849200	qPCR PVA/PVM set 192 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes/IC Mix_PVA/PVM/IC • Nuclease free water • Positive and negative controls 	192	
849610	qPCR PVA/PVM kit 96/10 includes: <ul style="list-style-type: none"> • qPCR PVA/PVM set 96 • Potato DNA/RNA rapid extraction set 10x for 96 extractions and a pool of 10 potatoes 	96	10
849625	qPCR PVA/PVM kit 96/25 includes: <ul style="list-style-type: none"> • qPCR PVA/PVM set 96 • Potato DNA/RNA rapid extraction set 25x for 96 extractions and a pool of 25 potatoes 	96	25
849210	qPCR PVA/PVM kit 192/10 includes: <ul style="list-style-type: none"> • qPCR PVA/PVM set 192 • Potato DNA/RNA rapid extraction set 10x for 192 extractions and a pool of 10 potatoes 	192	10

Art. No.	Product	No. of reactions	Amount of samples (potatoes)
849225	qPCR PVA/PVM kit 192/25 includes: <ul style="list-style-type: none"> • qPCR PVA/PVM set 192 • Potato DNA/RNA rapid extraction set 25x for 192 extractions and a pool of 25 potatoes 	192	25
879600	qPCR GRBV set 96 includes: <ul style="list-style-type: none"> • Taq Probe Master Mix • Primers/Probes Mix GRBV/IC • Nuclease free water 	96	
879200	qPCR GRBV set 192 includes: <ul style="list-style-type: none"> • Taq Probe Master Mix • Primers/Probes Mix GRBV/IC • Nuclease free water 	192	
889600*	qPCR ToBRFV set 96 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes Mix ToBRFV/COX • Nuclease free water 	96	
889200*	qPCR ToBRFV set 192 includes: <ul style="list-style-type: none"> • Taq Master Mix • RT Master Mix • Primers/Probes Mix ToBRFV/COX • Nuclease free water 	192	

* A ToBRFV RNA positive control (PC, Art. No. 880053, 30 µl) and a plant RNA negative control (NC, Art. No. 830043, 30 µl) is available separately upon request. It is mandatory to clarify before ordering, whether the ToBRFV PC may be imported into your country of destination.

ELISA Buffers



The ELISA buffers are supplied either as tablets or as concentrated liquids (5x or 10x concentrate). These formats are **easier** and **more accurate** to use than powder concentrates. There is no need for time-consuming dissolution of PVP and other chemicals or pH adjustment. Extraction, coating and conjugate buffers are available in different colours to avoid mistakes. ELISA buffers are available either individually or as part of complete kits (see page 17).

i

Additional Products

Product	Art. No.	Size	sufficient for
Coating buffer	110110	50 tablets	25000 assays
Coating buffer	110111	10 tablets	5000 assays
Coating buffer	110112	2 tablets	1000 assays
Conjugate buffer (10x)	110140	100 ml	5000 assays
Conjugate buffer (10x)	110142	20 ml	1000 assays
Substrate buffer* (5x)	110130	200 ml	5000 assays
Substrate buffer* (5x)	110132	40 ml	1000 assays
Washing buffer (PBST)	110150	10 tablets	10 x 1 liter
Washing buffer (PBST) powder	110152	50 g	5 liters
Washing buffer (PBST) powder	110153	250 g	25 liters
Washing buffer (PBST) powder	110154	1250 g	125 liters
Extraction buffer «General» (10x)	110120	500 ml	1000 samples**
Extraction buffer «Grapevine» (5x)	110123	2 x 500 ml	1000 samples**
Extraction buffer «Bulbs & Tubers» ***	110122	500 ml	1000 samples**
Extraction buffer «Caulimo» ****	110124	500 ml	1000 samples
Extraction buffer «Caulimo» ****	110224	250 ml	500 samples

* For use with pNPP (see below)

** Calculated with 5 ml extraction buffer per sample of 0.25 g

*** Composed of 500 ml extraction buffer «General» (10x) and 50 g egg albumin

**** Includes extraction buffer «Caulimo» (10x), urea, and egg albumin

Substrates

Product	Art. No.	Size	sufficient for
pNPP blister-packed tablets (20 mg)*	110180	10 tablets	1000 assays
pNPP blister-packed tablets (20 mg)*	110080	50 tablets	5000 assays
pNPP blister-packed tablets (5 mg)*	110185	10 tablets	250 assays
pNPP blister-packed tablets (5 mg)*	110085	50 tablets	1250 assays
pNPP bulk tablets (5 mg)**	50-80-01	100 tablets	2500 assays
pNPP bulk powder***	4165.8	1 x 5 g	25000 assays
pNPP bulk powder***	4165.9	1 x 25 g	125000 assays

* Single packed, providing optimal protection and long-term stability (shelf life: 3 years at -20°C)


** Bulk reagent (shelf life: 3 years at -20°C)

*** Excellent quality bulk substrate; needs a balance with 0.1 mg accuracy (shelf life: 5 years at -20°C)

Chemicals

Product	Art. No.	Size
BSA standard grade	700110	75 g
Egg albumin	5461.250	250 g
Egg albumin	5461.1000	1 kg

Consumables



ELISA Plates

Product	Art. No.	Size
Nunc-Immuno Plates MaxiSorp F96, certified	520000	10
Standard packaging: 6 x 10 plates	520000	60
ELISA plate covers (for multiple use)*	520001	10
Rubber sheets for ELISA plates (for multiple use)*	500003	10

*ELISA plates are covered during incubation with a sheet of parafilm which is pressed down firmly by a rubber sheet inside a plate cover. The plate cover is fixed to the plate with rubber bands. This prevents evaporation during incubation.



Reservoirs

Product	Art. No.	Size
Single reservoir 70 ml (polypropylene)	330.01	20
8-channel reservoir 40 ml (polystyrene)	330.08.9	10
12-channel reservoir 40 ml (polystyrene)	330.12.9	10

BIOREBA Testing Service Laboratory for Plant Pathogen Diagnostics



Association of Swiss Laboratories
Verband Schweizer Laboratorien
Association des Laboratoires Suisses
Associazione dei Laboratori Svizzeri


BIOREBA Testing Service offers you its expertise in diagnostic assays for the detection of plant pathogens (viruses, bacteria, fungi and phytoplasmas).

More than 30 years of experience with analysis of samples of most agricultural sectors and employment of different state-of-the-art methods and technologies guarantee the **highest** possible accuracy of test results.

Each case is handled in **close** communication with the client to determine the optimal sampling and transport of samples, as well as the most appropriate testing method and conditions. This ensures that each customer receives the right service for their needs in the most cost effective and timely manner.

Our testing service complies with the strict rules and quality standards of **ISO/IEC 17025** and confidentiality is maintained in all communications and data processing.

For further information please visit our website www.bioreba.com or contact us by telephone or email at admin@bioreba.ch.



Prices: Please contact us now for a quote. Our pricing is based on sample volume, sample types, number of pathogens to be tested per sample, and testing method.



Quality Management

Our **vision** is to provide products and services of outstanding quality to our customers around the world. Our Quality Management System is an integral part of all our operations. Our value added processes are constantly being improved. Within the framework of our quality management system, we are constantly looking for ways to further develop our company.

BIOREBA has been ISO certified for many years. Our quality management system in accordance with **ISO 9001:2015** applies not only to the production, research and development stages, but also to the worldwide distribution of our products and the testing service for customers.

In addition, we share our own experience with our customers in our **ISO/IEC 17025** accredited testing laboratory and support them in process validation or other quality-related issues. BIOREBA is an official testing laboratory for phytodiagnostics. As such, we participate in numerous Test Performance Studies (TPS).



THE INTERNATIONAL CERTIFICATION NETWORK

CERTIFICATE

SQS has issued an IQNet recognized certificate that the organization:

BIOREBA AG
 Christoph Merian-Ring 7
 4153 Reinach
 Switzerland

has implemented and maintains a
Management System
for the following scope:

Agro-Diagnostics: Research and development, production and worldwide distribution of tests, test systems and equipment, as well as testing service for customers

which fulfills the requirements of the following standard(s):

ISO 9001:2015

Issued on: 2021-02-16
Expires on: 2024-02-15

This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document

Registration Number: CH-35633





Alex Stoichitou
President of IQNet

Felix Müller
CEO SQS

IQNet Partners:
 AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISO Italy
 CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany EAGLE Certification Group USA
 FCV Brazil FONDAZIONE Veneranda ICA-INTZ Colombia Inspecta Sertifiointi Oy Finland INTECO Costa Rica
 IRAM Argentina JQA Japan KPO Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland
 NYCE-SICE Mexico PCBC Poland Quality Austria Austria RS Russia SII Israel SIQ Slovenia
 SIRIM QAS International Malaysia SQS Switzerland SKAC Romania TEST StPetersburg Russia TSE Turkey YUQS Serbia

* The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under www.iqnet-certification.com



Schweizerische Eidgenossenschaft
Confédération suisse
Confederazione Svizzera
Confederaziun svizra

Swiss Confederation

Federal Department of Economic Affairs,
Education and Research EAER
State Secretariat for Economic Affairs SECO
Swiss Accreditation Service SAS

Based on the Accreditation and Designation Ordinance dated 17 June 1996 and on the advice of the Federal Accreditation Commission, the Swiss Accreditation Service (SAS) grants to

BIOREBA AG
 Test service
 Plant pathogen diagnostics
 Chr. Merian-Ring 7
 4153 Reinach



Period of accreditation:
28.05.2020 until 27.05.2025
 (1st accreditation: 28.05.2010)

the accreditation as

Testing laboratory for serological and molecular biological analyses of plant pathogens

International standard: ISO/IEC 17025:2017
 Swiss standard: SN EN ISO/IEC 17025:2018

3003 Berne, 12.02.2020
 Swiss Accreditation Service SAS



Head of SAS
 Konrad Flück

SAS is a signatory of the multilateral agreements of the European co-operation for Accreditation (EA) for the fields of testing, calibration, inspection and certification of management systems, certification of personnel and certification of products, processes and services, of the International Accreditation Forum (IAF) for the fields of certification of management systems and certification of products, processes and services and of the International Laboratory Accreditation Cooperation (ILAC) for the fields of testing and calibration.

Worldwide distribution



BIOREBA is committed to providing our sample preparation products, plant pathogen diagnostic tests and ancillary products on a global scale. We supply our products worldwide. In addition, we have established long-term strategic relationships with local partners in **more than 50 countries**.

You can find a list of our distributors on our website:

www.bioreba.com

For direct sales support, please contact your local distributor or BIOREBA at our headquarters in Switzerland.



Mission Statement

“With the detection of plant diseases, we enable to promote a stable, consistent food supply for humans and animals **worldwide.**”



*Your Partner in
Agro-Diagnostics*
www.bioreba.com

Version: 02.2023 - 10.07.2023
Adaptions from last version: CLRv-e deleted.

