

Product Information

Substrates



Para-nitrophenyl-phosphate (pNPP)	Art. No.	Size	sufficient for
pNPP blister-packed tablets (20 mg)*	110180	10 tablets	1000 assays
pNPP blister-packed tablets (20 mg)*	110080	50 tablets	5000 assays
pNPP blister-packed tablets (5 mg)*	110185	10 tablets	250 assays
pNPP blister-packed tablets (5 mg)*	110085	50 tablets	1250 assays
pNPP bulk tablets (5 mg)**	50-80-01	100 tablets	2500 assays
pNPP bulk powder***	4165.8	1 x 5 g	25000 assays
pNPP bulk powder***	4165.9	1 x 25 g	125000 assays

* Single packed, providing optimal protection and long-term stability (shelf life: 3 years at -20°C)

** Bulk reagent with limited shelf life (shelf life: 1 year, at 4-6°C)

*** Excellent quality bulk substrate; needs a balance with 0.1 mg accuracy (shelf life: 2 years at 4-6°C)